






# GYNAECOLOGY



GYNAECOLOGY

# ABOUT US



SINCE 20 YEARS



EXPORT IN  
+70 COUNTRIES



24H/SERVICE



+ 10.000  
PRODUCTS



LOCATED IN THE  
WORLD CENTER  
FOR MEDICAL  
TECHNOLOGY

## CERTIFICATED

STEMA Medizintechnik GmbH is certified in accordance with the highest standards:

- EN ISO 13485:2016 + AC:2018 + A11:2021- ISO 13485:2016
- EU 2017/745 (**MDR**)
- Council Directive 93/42/EEC Annex II excluding Section 4
- Council Directive 93/42/EEC Annex V excluding Section 3
- FDA

These standards create an essential base for our company and our partners. They also symbolize the commitment and strive to keep up the quality of our products.



# PRODUCT RANGE

The current product range of STEMA Medizintechnik GmbH includes more than 10.000 different items.

We cover the following areas in surgery with high quality products:

ARTHROSCOPY  
BRONCHOSCOPY  
CARDIAC MIS  
CONTAINER-SYSTEMS  
GASTROENTEROLOGY  
GYNAECOLOGY  
HF-SURGERY  
LAMINECTOMY  
LAPAROSCOPY  
MICRO MOTOR SYSTEMS  
PROCTOSCOPY  
SPINE-SURGERY  
THORACOSCOPY  
TITANIUM PLATING SYSTEM  
UNITS  
UROLOGY



# CONTENT

## GYNAECOLOGY

### PRODUCT

### PAGE

#### 1. HYSTEROSCOPY 4.0MM

|                             |    |
|-----------------------------|----|
| 1.1 Cystoscopes 4.0mm       | 8  |
| 1.2 Diagnostic Hysteroscopy | 9  |
| 1.3 Operating Hysteroscopy  | 11 |

#### 2. SLIM HYSTEROSCOPY 2.9MM

|                                 |    |
|---------------------------------|----|
| 2.1 Cystoscopes 2.9mm           | 17 |
| 2.2 Diagnostic Hysteroscopy     | 18 |
| 2.3 Operating Hysteroscopy      | 20 |
| 2.4 High-Flow Operating sheaths | 23 |

#### 3. MINI HYSTEROSCOPY 2.0MM

|                                 |    |
|---------------------------------|----|
| 3.1 Cystoscopes 2.0mm           | 26 |
| 3.2 Diagnostic Hysteroscopy     | 27 |
| 3.3 Operating Hysteroscopy      | 29 |
| 3.4 High-Flow Operating sheaths | 30 |

#### 4. INSTRUMENTS & ELECTRODES

|                                    |    |
|------------------------------------|----|
| 4.1 Forceps, scissors, instruments | 33 |
| 4.2 Electrodes & connecting cables | 40 |
| 4.3 Accessories                    | 45 |

#### 5. SALINE RESECTOSCOPY

|                              |    |
|------------------------------|----|
| 5.1 Saline Resectoscopy      | 59 |
| 5.2 Saline Slim Resectoscopy | 65 |
| 5.3 Saline Mini Resectoscopy | 70 |
| 5.4 Connecting Cables        | 75 |

# CONTENT

## GYNAECOLOGY

### PRODUCT

### PAGE

#### 6. HYSTEROSCOPY SETS

|                                                       |    |
|-------------------------------------------------------|----|
| 6.1 Diagnostic Hysteroscopy set                       | 77 |
| 6.2 Operating Hysteroscopy set                        | 78 |
| 6.3 Out-Patient Atraumatic Operating Hysteroscopy set | 79 |
| 6.4 Saline Resectoscopy set                           | 80 |

#### 7. INDEX

|           |    |
|-----------|----|
| 7.1 Index | 83 |
|-----------|----|

# 1. HYSTEROSCOPY 4.0MM

1.1 | CYSTOSCOPES 4.0MM


1.2 | DIAGNOSTIC HYSTEROSCOPY

1.3 | OPERATING HYSTEROSCOPY

# 1. HYSTEROSCOPY

## COMPLETE DIAGNOSTIC HYSTEROSCOPE

4.0MM




### 1 CYSTOSCOPES



### 2 SHEATHS



### 3 OBTURATORS



## 1.1 | HYSTEROSCOPY

### Ø 4.0MM CYSTOSCOPES AUTOCLAVABLE



| Viewing angle                                                                           | Endoscope                  | Diameter | Working length | Ordering no. |
|-----------------------------------------------------------------------------------------|----------------------------|----------|----------------|--------------|
| <br>0° | Cystoscope autoclavable    | Ø 4.0mm  | 302mm          | AA-002-001   |
|                                                                                         | HD Cystoscope autoclavable | Ø 4.0mm  | 302mm          | AA-002-011   |

| Viewing angle                                                                              | Endoscope                  | Diameter | Working length | Ordering no. |
|--------------------------------------------------------------------------------------------|----------------------------|----------|----------------|--------------|
| <br>12° | Cystoscope autoclavable    | Ø 4.0mm  | 302mm          | AA-002-002   |
|                                                                                            | HD Cystoscope autoclavable | Ø 4.0mm  | 302mm          | AA-002-012   |


| Viewing angle                                                                              | Endoscope                  | Diameter | Working length | Ordering no. |
|--------------------------------------------------------------------------------------------|----------------------------|----------|----------------|--------------|
| <br>30° | Cystoscope autoclavable    | Ø 4.0mm  | 302mm          | AA-002-003   |
|                                                                                            | HD Cystoscope autoclavable | Ø 4.0mm  | 302mm          | AA-002-013   |

| Viewing angle                                                                              | Endoscope                  | Diameter | Working length | Ordering no. |
|--------------------------------------------------------------------------------------------|----------------------------|----------|----------------|--------------|
| <br>70° | HD Cystoscope autoclavable | Ø 4.0mm  | 302mm          | AA-002-014   |



## 1.2 | DIAGNOSTIC HYSTEROSCOPY

### SHEATHS & OBTURATOR



| Ordering no.      | Diameter | WL    | Description                                                      |
|-------------------|----------|-------|------------------------------------------------------------------|
| <b>AN-003-001</b> | Ø 5.1mm  | 277mm | <b>Diagnostic</b> hysteroscope sheath for Ø 4.0mm, 0° endoscope  |
| <b>AN-003-002</b> | Ø 5.1mm  | 277mm | <b>Diagnostic</b> hysteroscope sheath for Ø 4.0mm, 30° endoscope |



| Ordering no.      | Description                                                               |
|-------------------|---------------------------------------------------------------------------|
| <b>AN-003-003</b> | Obturator for <b>diagnostic</b> hysteroscope sheath for Ø 4.0mm endoscope |

## 1.2 | CONTINUOUS-FLOW DIAGNOSTIC HYSTEROSCOPY SHEATHS & OBTURATOR



| Ordering no.      | Diameter | WL    | Description                                                                            |
|-------------------|----------|-------|----------------------------------------------------------------------------------------|
| <b>AN-003-103</b> | Ø 6.5mm  | 246mm | <b>Outer-Diagnostic-Continuous-Flow</b> hysteroscope sheath for <b>4.0mm</b> endoscope |




| Ordering no.      | Diameter | WL    | Description                                                                                |
|-------------------|----------|-------|--------------------------------------------------------------------------------------------|
| <b>AN-003-101</b> | Ø 5.4mm  | 261mm | <b>Inner-Diagnostic-Continuous-Flow</b> hysteroscope sheath for <b>4.0mm/0°</b> endoscope  |
| <b>AN-003-102</b> | Ø 5.4mm  | 261mm | <b>Inner-Diagnostic-Continuous-Flow</b> hysteroscope sheath for <b>4.0mm/30°</b> endoscope |



| Ordering no.      | Description                                                                                          |
|-------------------|------------------------------------------------------------------------------------------------------|
| <b>AN-003-104</b> | Obturator for <b>Outer-Diagnostic-Continuous-Flow</b> hysteroscope sheath for <b>4.0mm</b> endoscope |

## 1.3 | CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS



| Ordering no.      | Diameter | WL    | Description                                                                                                                             |
|-------------------|----------|-------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>AN-003-201</b> | Ø 8.0mm  | 193mm | <b>Standard OP-Continuous-Flow</b> hysteroscope sheath for <b>0°/12°</b> endoscope and <b>7 Fr. flexible and semi-rigid instruments</b> |
| <b>AN-003-202</b> | Ø 8.0mm  | 193mm | <b>Standard OP-Continuous-Flow</b> hysteroscope sheath for <b>30°</b> endoscope and <b>7 Fr. flexible and semi-rigid instruments</b>    |

### 1.3 | CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS & OBTURATOR



| Ordering no.      | Diameter | WL    | Description                                                                                                                          |
|-------------------|----------|-------|--------------------------------------------------------------------------------------------------------------------------------------|
| <b>AN-003-302</b> | Ø 8.0mm  | 180mm | <b>Outer-OP-Continuous-Flow</b> hysteroscope sheath for <b>12°/30</b> endoscope and <b>7 Fr. flexible and semi-rigid instruments</b> |



| Ordering no.      | Diameter | WL    | Description                                                                                                                          |
|-------------------|----------|-------|--------------------------------------------------------------------------------------------------------------------------------------|
| <b>AN-003-301</b> | Ø 8.0mm  | 180mm | <b>Inner-OP-Continuous-Flow</b> hysteroscope sheath for <b>12°/30</b> endoscope and <b>7 Fr. flexible and semi-rigid instruments</b> |



| Ordering no.      | Diameter | WL    | Description                                                       |
|-------------------|----------|-------|-------------------------------------------------------------------|
| <b>AN-003-303</b> | Ø 8.0mm  | 180mm | Obturator for <b>Outer-OP-Continuous-Flow</b> hysteroscope sheath |

### 1.3 | CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS



| Ordering no.      | Diameter | WL    | Description                                                                                                                   |
|-------------------|----------|-------|-------------------------------------------------------------------------------------------------------------------------------|
| <b>AN-005-302</b> | Ø 6.5mm  | 192mm | <b>Outer-OP-Continuous-Flow</b> hysteroscope sheath for <b>30°</b> endoscope and <b>5 Fr. flexible/semi-rigid instruments</b> |
| <b>AN-005-304</b> | Ø 7.0mm  | 192mm | <b>Outer-OP-Continuous-Flow</b> hysteroscope sheath for <b>30°</b> endoscope and <b>7 Fr. flexible/semi-rigid instruments</b> |



| Ordering no.      | Diameter | WL    | Description                                                                                                                   |
|-------------------|----------|-------|-------------------------------------------------------------------------------------------------------------------------------|
| <b>AN-005-301</b> | Ø 6.5mm  | 192mm | <b>Inner-OP-Continuous-Flow</b> hysteroscope sheath for <b>30°</b> endoscope and <b>5 Fr. flexible/semi-rigid instruments</b> |
| <b>AN-005-303</b> | Ø 7.0mm  | 192mm | <b>Inner-OP-Continuous-Flow</b> hysteroscope sheath for <b>30°</b> endoscope and <b>7 Fr. flexible/semi-rigid instruments</b> |

### 1.3 | CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS & OBTURATOR



| Ordering no.      | Diameter | WL    | Description                                                                                                    |
|-------------------|----------|-------|----------------------------------------------------------------------------------------------------------------|
| <b>AN-004-401</b> | Ø 8.0mm  | 180mm | <b>Outer-OP-Continuous-Flow</b> hysteroscope sheath for <b>7 Fr. flexible and semi-rigid 12°/30° endoscope</b> |



| Ordering no.      | Diameter | WL    | Description                                    |
|-------------------|----------|-------|------------------------------------------------|
| <b>AN-004-402</b> | Ø 8.0mm  | 180mm | Working Element with <b>Albarran Deflector</b> |



| Ordering no.      | Diameter | WL    | Description                                                                                  |
|-------------------|----------|-------|----------------------------------------------------------------------------------------------|
| <b>AN-004-403</b> | Ø 8.0mm  | 180mm | Obturator for <b>Outer-OP-Continuous-Flow</b> hysteroscope sheath for <b>4.0mm endoscope</b> |

## 2. SLIM HYSTEROSCOPY 2.9MM

2.1 | CYSTOSCOPES 2.9MM

2.2 | DIAGNOSTIC HYSTEROSCOPY

2.3 | OPERATING HYSTEROSCOPY

2.4 | HIGH-FLOW OPERATING SHEATHS

## 2. SLIM HYSTEROSCOPY 2.9MM


2.9MM



### 1 CYSTOSCOPES



### 2 SHEATHS



### 3 OBTURATORS





## 2.1 | SLIM HYSTEROSCOPY

### Ø 2.9MM CYSTOSCOPES AUTOCLAVABLE



| Viewing angle                                                                     | Endoscope                         | Diameter | Working length | Ordering no.      |
|-----------------------------------------------------------------------------------|-----------------------------------|----------|----------------|-------------------|
|  | Cystoscope autoclavable           | Ø 2.9mm  | 302mm          | <b>AA-003-001</b> |
|                                                                                   | <b>HD</b> Cystoscope autoclavable | Ø 2.9mm  | 302mm          | <b>AA-003-011</b> |
| <b>0°</b>                                                                         |                                   |          |                |                   |

| Viewing angle                                                                       | Endoscope                         | Diameter | Working length | Ordering no.      |
|-------------------------------------------------------------------------------------|-----------------------------------|----------|----------------|-------------------|
|  | <b>HD</b> Cystoscope autoclavable | Ø 2.9mm  | 302mm          | <b>AA-003-014</b> |
|                                                                                     |                                   |          |                |                   |
| <b>12°</b>                                                                          |                                   |          |                |                   |

| Viewing angle                                                                       | Endoscope                         | Diameter | Working length | Ordering no.      |
|-------------------------------------------------------------------------------------|-----------------------------------|----------|----------------|-------------------|
|  | Cystoscope autoclavable           | Ø 2.9mm  | 302mm          | <b>AA-003-002</b> |
|                                                                                     | <b>HD</b> Cystoscope autoclavable | Ø 2.9mm  | 302mm          | <b>AA-003-012</b> |
| <b>30°</b>                                                                          |                                   |          |                |                   |

## 2.2 | SLIM DIAGNOSTIC HYSTEROSCOPY SHEATHS & OBTURATOR



| Ordering no.      | Diameter | WL    | Description                                                      |
|-------------------|----------|-------|------------------------------------------------------------------|
| <b>AN-002-001</b> | Ø 4.0mm  | 277mm | <b>Diagnostic</b> hysteroscope sheath for Ø 2.9mm, 0° endoscope  |
| <b>AN-002-002</b> | Ø 4.0mm  | 277mm | <b>Diagnostic</b> hysteroscope sheath for Ø 2.9mm, 30° endoscope |



| Ordering no.      | Description                                                             |
|-------------------|-------------------------------------------------------------------------|
| <b>AN-002-003</b> | Obturator for <b>diagnostic</b> hysteroscope sheath for 2.9mm endoscope |

## 2.2 | SLIM CONTINUOUS-FLOW DIAGNOSTIC HYSTEROSCOPY SHEATHS & OBTURATOR



| Ordering no.      | Diameter | WL    | Colour | Description                                                                            |
|-------------------|----------|-------|--------|----------------------------------------------------------------------------------------|
| <b>AN-002-103</b> | Ø 5mm    | 246mm | White  | <b>Outer-Diagnostic-Continuous-Flow</b> hysteroscope sheath for <b>2.9mm</b> endoscope |



| Ordering no.      | Diameter | WL    | Colour | Description                                                                                |
|-------------------|----------|-------|--------|--------------------------------------------------------------------------------------------|
| <b>AN-002-101</b> | Ø 3.7mm  | 261mm | White  | <b>Inner-Diagnostic-Continuous-Flow</b> hysteroscope sheath for <b>2.9mm/0°</b> endoscope  |
| <b>AN-002-102</b> | Ø 3.7mm  | 261mm | White  | <b>Inner-Diagnostic-Continuous-Flow</b> hysteroscope sheath for <b>2.9mm/30°</b> endoscope |



| Ordering no.      | Description                                                                                          |
|-------------------|------------------------------------------------------------------------------------------------------|
| <b>AN-002-104</b> | Obturator for <b>Outer-Diagnostic-Continuous-Flow</b> hysteroscope sheath for <b>2.9mm</b> endoscope |

## 2.3 | SLIM CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS



| Ordering no.      | Diameter | WL    | Description                                                                                                                            |
|-------------------|----------|-------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>AN-002-201</b> | Ø 6.0mm  | 193mm | <b>Standard OP-Continuous-Flow</b> hysteroscope sheath<br><b>0°/12° endoscope</b> and <b>5 Fr. flexible and semi-rigid instruments</b> |
| <b>AN-002-202</b> | Ø 6.0mm  | 193mm | <b>Standard OP-Continuous-Flow</b> hysteroscope sheath<br><b>30° endoscope</b> and <b>7 Fr. flexible and semi-rigid instruments</b>    |

## 2.3 | SLIM CONTINOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS & OBTURATOR



| Ordering no.      | Diameter | WL    | Description                                                                                                             |
|-------------------|----------|-------|-------------------------------------------------------------------------------------------------------------------------|
| <b>AN-002-302</b> | Ø 5.0mm  | 192mm | <b>Outer-OP-Continuous-Flow</b> hysteroscope sheath for 30° endoscope and 5 Fr. flexible/semi-rigid instruments, white  |
| <b>AN-002-308</b> | Ø 6.3mm  | 192mm | <b>Outer-OP-Continuous-Flow</b> hysteroscope sheath for 30° endoscope and 9 Fr. flexible/semi-rigid instruments, yellow |



| Ordering no.      | Diameter | WL    | Description                                                                                                             |
|-------------------|----------|-------|-------------------------------------------------------------------------------------------------------------------------|
| <b>AN-002-301</b> | Ø 5.0mm  | 192mm | <b>Inner-OP-Continuous-Flow</b> hysteroscope sheath for 30° endoscope and 5 Fr. flexible/semi-rigid instruments, white  |
| <b>AN-002-307</b> | Ø 6.3mm  | 192mm | <b>Inner-OP-Continuous-Flow</b> hysteroscope sheath for 30° endoscope and 9 Fr. flexible/semi-rigid instruments, yellow |

## 2.3 | SLIM CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS



| Ordering no.      | Diameter       | WL    | Description                                                                                                            |
|-------------------|----------------|-------|------------------------------------------------------------------------------------------------------------------------|
| <b>AN-002-312</b> | <b>Ø 5.0mm</b> | 216mm | <b>Outer-OP-Continuous-Flow</b> hysteroscope sheath for <b>30°</b> endoscope and 5 Fr. flexible/semi-rigid instruments |



| Ordering no.      | Diameter       | WL    | Description                                                                                                            |
|-------------------|----------------|-------|------------------------------------------------------------------------------------------------------------------------|
| <b>AN-002-313</b> | <b>Ø 5.0mm</b> | 216mm | <b>Inner-OP-Continuous-Flow</b> hysteroscope sheath for <b>30°</b> endoscope and 5 Fr. flexible/semi-rigid instruments |

## 2.4 | SLIM HIGH-FLOW CONTINUOUS-FLOW OPERATING SHEATHS



| Ordering no.      | Diameter | Description                                                                                                                                                                       |
|-------------------|----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>AN-002-310</b> | Ø 5.5mm  | <b>High-Flow Continuous-Flow Operating</b> hysteroscope sheath for <b>5 Fr. flexible and semi-rigid instruments for use with 2.9mm endoscope (12°/30°)</b> , single sheath system |



| Ordering no.      | Diameter | Description                                                                                                                                                                                             |
|-------------------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>AN-002-311</b> | Ø 5.5mm  | <b>High-Flow Continuous-Flow Operating</b> hysteroscope sheath, for <b>5 Fr. flexible and semi-rigid instruments for use with 2.9mm endoscope (12°/30°)</b> , single sheath system <b>with stopcock</b> |

## 3. MINI HYSTEROSCOPY 2.0MM

**3.1 | CYSTOSCOPES 2.0MM**

**3.2 | DIAGNOSTIC SHEATHS & OBTURATORS**


**3.3 | OPERATING SHEATHS & OBTURATORS**

**3.4 | HIGH-FLOW OPERATING SHEATHS**



### 3. MINI HYSTEROSCOPY 2.0MM

2MM



#### 1 CYTOSCOPES



#### 2 SHEATHS



#### 3 OBTURATORS



### 3.1 | MINI HYSTEROSCOPY

#### Ø 2.0MM CYSTOSCOPES AUTOCLAVABLE



| Viewing angle                                                                           | Endoscope               | Diameter | Working length | Ordering no. |
|-----------------------------------------------------------------------------------------|-------------------------|----------|----------------|--------------|
| <br>0° | Cystoscope autoclavable | Ø 2.0mm  | 260mm          | AA-004-001   |

| Viewing angle                                                                              | Endoscope               | Diameter | Working length | Ordering no. |
|--------------------------------------------------------------------------------------------|-------------------------|----------|----------------|--------------|
| <br>12° | Cystoscope autoclavable | Ø 2.0mm  | 260mm          | AA-004-003   |

| Viewing angle                                                                              | Endoscope               | Diameter | Working length | Ordering no. |
|--------------------------------------------------------------------------------------------|-------------------------|----------|----------------|--------------|
| <br>30° | Cystoscope autoclavable | Ø 2.0mm  | 260mm          | AA-004-002   |

### 3.2 | MINI DIAGNOSTIC HYSTEROSCOPY SHEATHS & OBTURATOR



| Ordering no.      | Diameter | WL    | Description                                                      |
|-------------------|----------|-------|------------------------------------------------------------------|
| <b>AN-001-001</b> | Ø 2.5mm  | 277mm | <b>Diagnostic</b> hysteroscope sheath for Ø 2.0mm, 0° endoscope  |
| <b>AN-001-002</b> | Ø 2.5mm  | 277mm | <b>Diagnostic</b> hysteroscope sheath for Ø 2.0mm, 30° endoscope |



| Ordering no.      | Description                                                             |
|-------------------|-------------------------------------------------------------------------|
| <b>AN-001-003</b> | Obturator for <b>diagnostic</b> hysteroscope sheath for 2.0mm endoscope |

### 3.2 | MINI CONTINUOUS-FLOW DIAGNOSTIC HYSTEROSCOPY SHEATHS & OBTURATOR



| Ordering no.      | Diameter | WL    | Colour | Description                                                                              |
|-------------------|----------|-------|--------|------------------------------------------------------------------------------------------|
| <b>AN-001-103</b> | Ø 3.4mm  | 204mm | Red    | <b>Outer-Diagnostic-Continuous-Flow</b> hysteroscope sheath for <b>Ø 2.0mm endoscope</b> |



| Ordering no.      | Diameter | WL    | Colour | Description                                                                                   |
|-------------------|----------|-------|--------|-----------------------------------------------------------------------------------------------|
| <b>AN-001-101</b> | Ø 2.6mm  | 219mm | Red    | <b>Inner-Diagnostic-Continuous-Flow</b> hysteroscope sheath for <b>Ø 2.0mm, 0° endoscope</b>  |
| <b>AN-001-102</b> | Ø 2.6mm  | 219mm | Red    | <b>Inner-Diagnostic-Continuous-Flow</b> hysteroscope sheath for <b>Ø 2.0mm, 30° endoscope</b> |



| Ordering no.      | Description                                                                                          |
|-------------------|------------------------------------------------------------------------------------------------------|
| <b>AN-001-104</b> | Obturator for <b>Outer-Diagnostic-Continuous-Flow</b> hysteroscope sheath for <b>2.0mm endoscope</b> |

### 3.3 | MINI CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SHEATHS




| Ordering no.      | Diameter | WL    | Description                                                                                                                             |
|-------------------|----------|-------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>AN-001-302</b> | Ø 4.3mm  | 181mm | <b>Outer-OP-Continuous-Flow</b> hysteroscope sheath for <b>2.0mm/30°</b> endoscope and <b>5 Fr. flexible and semi-rigid instruments</b> |




| Ordering no.      | Diameter | WL    | Description                                                                                                                                     |
|-------------------|----------|-------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>AN-001-301</b> | Ø 4.3mm  | 181mm | <b>Inner-Diagnostic-Continuous-Flow</b> hysteroscope sheath for <b>2.0mm/30°</b> endoscope and <b>5 Fr. flexible and semi-rigid instruments</b> |

### 3.4 | HIGH-FLOW CONTINUOUS-FLOW OPERATING SHEATHS



| Ordering no.      | Diameter | Description                                                                                                                                                                       |
|-------------------|----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>AN-001-310</b> | Ø 4.3mm  | <b>High-Flow Continuous-Flow Operating</b> hysteroscope sheath for <b>5 Fr. flexible and semi-rigid instruments for use with 2.0mm endoscope (12°/30°)</b> , single sheath system |



| Ordering no.      | Diameter | Description                                                                                                                                                                                             |
|-------------------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>AN-001-311</b> | Ø 4.3mm  | <b>High-Flow Continuous-Flow Operating</b> hysteroscope sheath, for <b>5 Fr. flexible and semi-rigid instruments for use with 2.0mm endoscope (12°/30°)</b> , single sheath system <b>with stopcock</b> |

## 4. INSTRUMENTS & ACCESSORIES

4.1 | FORCEPS, SCISSORS, INSTRUMENTS

4.2 | ELECTRODES & CONNECTING CABLES

4.3 | ACCESSORIES & HYSTEROSCOPY INSTRUMENTS

# 4. INSTRUMENTS & ACCESSORIES


## 1 INSTRUMENTS



## 2 ELECTRODES



## 3 INTRAUTERINE PROBE




## 4 RETRACTORS





## 4.1 | FORCEPS, SCISSORS, INSTRUMENTS


### FORCEPS FLEXIBLE | SEMI-RIGID | RIGID




|                                                                                    | Description                                                   | WL   | Diameter | French | Flexible   | Semi-rigid | Rigid        |
|------------------------------------------------------------------------------------|---------------------------------------------------------------|------|----------|--------|------------|------------|--------------|
|    | Biopsy forceps, single action                                 | 40cm | Ø 2.3mm  | 7 Fr.  | AC-010-004 | AC-010-204 |              |
|   | Grasping forceps, single action                               | 40cm | Ø 2.3mm  | 7 Fr.  | AC-011-004 | AC-011-204 |              |
|  | Biopsy forceps oval, double action                            | 40cm | Ø 1.0mm  | 3 Fr.  | AC-001-001 | AC-001-201 | AC-001-001-R |
|                                                                                    |                                                               | 40cm | Ø 1.2mm  | 4 Fr.  | AC-001-008 | AC-001-208 | AC-001-008-R |
|                                                                                    |                                                               | 40cm | Ø 1.6mm  | 5 Fr.  | AC-001-002 | AC-001-202 | AC-001-002-R |
|                                                                                    |                                                               | 40cm | Ø 2.1mm  | 6 Fr.  | AC-001-003 | AC-001-203 | AC-001-003-R |
|                                                                                    |                                                               | 40cm | Ø 2.3mm  | 7 Fr.  | AC-001-004 | AC-001-204 | AC-001-004-R |
|                                                                                    |                                                               | 40cm | Ø 3.0mm  | 9 Fr.  | AC-001-005 | AC-001-205 | AC-001-005-R |
|  | Biopsy forceps oval, serrated, double action                  | 40cm | Ø 1.0mm  | 3 Fr.  | AC-002-001 | AC-002-201 | AC-002-001-R |
|                                                                                    |                                                               | 40cm | Ø 1.6mm  | 5 Fr.  | AC-002-002 | AC-002-202 | AC-002-002-R |
|                                                                                    |                                                               | 40cm | Ø 2.1mm  | 6 Fr.  | AC-002-003 | AC-002-203 | AC-002-003-R |
|                                                                                    |                                                               | 40cm | Ø 2.3mm  | 7 Fr.  | AC-002-004 | AC-002-204 | AC-002-004-R |
|                                                                                    |                                                               | 40cm | Ø 3.0mm  | 9 Fr.  | AC-002-005 | AC-002-205 | AC-002-005-R |
|  | Biopsy forceps oval, serrated, teeth interlock, double action | 40cm | Ø 1.6mm  | 5 Fr.  | AC-015-002 | AC-015-202 | AC-015-002-R |
|                                                                                    |                                                               | 40cm | Ø 2.1mm  | 6 Fr.  |            | AC-015-203 | AC-015-003-R |
|                                                                                    |                                                               | 40cm | Ø 2.3mm  | 7 Fr.  | AC-015-004 | AC-015-204 | AC-015-004-R |
|  | Grasping forceps, double action                               | 40cm | Ø 1.2mm  | 4 Fr.  | AC-017-001 | AC-017-201 | AC-017-001-R |
|                                                                                    |                                                               | 40cm | Ø 2.3mm  | 7 Fr.  | AC-017-004 | AC-017-204 | AC-017-004-R |
|  | Grasping forceps, alligator jaw, double action                | 40cm | Ø 1.0mm  | 3 Fr.  | AC-006-001 | AC-006-201 | AC-006-001-R |
|                                                                                    |                                                               | 40cm | Ø 1.2mm  | 4 Fr.  | AC-006-003 | AC-006-203 | AC-006-003-R |
|                                                                                    |                                                               | 40cm | Ø 1.6mm  | 5 Fr.  | AC-006-002 | AC-006-202 | AC-006-002-R |
|                                                                                    |                                                               | 40cm | Ø 2.3mm  | 7 Fr.  | AC-006-004 | AC-006-204 | AC-006-004-R |
|                                                                                    |                                                               | 40cm | Ø 3.0mm  | 9 Fr.  | AC-006-005 | AC-006-205 | AC-006-005-R |

## 4.1 | FORCEPS, SCISSORS, INSTRUMENTS


### FORCEPS FLEXIBLE | SEMI-RIGID | RIGID



|                                                                                    | Description                                      | WL   | Diameter | French | Flexible          | Semi-rigid                                                                            | Rigid               |
|------------------------------------------------------------------------------------|--------------------------------------------------|------|----------|--------|-------------------|---------------------------------------------------------------------------------------|---------------------|
|    | Grasping forceps, DJ, double action              | 40cm | Ø 2.3mm  | 7 Fr.  | <b>AC-008-004</b> | <b>AC-008-204</b>                                                                     | <b>AC-008-004-R</b> |
|                                                                                    |                                                  | 40cm | Ø 3.0mm  | 9 Fr.  | <b>AC-008-005</b> | <b>AC-008-205</b>                                                                     | <b>AC-008-005-R</b> |
|   | Grasping forceps teeth double action             | 40cm | Ø 1.2mm  | 4 Fr.  | <b>AC-016-001</b> | <b>AC-016-201</b>                                                                     | <b>AC-016-001-R</b> |
|                                                                                    |                                                  | 40cm | Ø 1.6mm  | 5 Fr.  | <b>AC-016-002</b> | <b>AC-016-202</b>                                                                     | <b>AC-016-002-R</b> |
|  | Grasping forceps teeth interlock, double action  | 40cm | Ø 1.6mm  | 5 Fr.  | <b>AC-020-002</b> | <b>AC-020-202</b>                                                                     | <b>AC-020-002-R</b> |
|                                                                                    |                                                  | 40cm | Ø 2.3mm  | 7 Fr.  | <b>AC-020-004</b> | <b>AC-020-204</b>                                                                     | <b>AC-020-004-R</b> |
|  | Grasping forceps rat tooth, double action        | 40cm | Ø 1.0mm  | 3 Fr.  | <b>AC-019-001</b> | <b>AC-019-201</b>                                                                     | <b>AC-019-001-R</b> |
|                                                                                    |                                                  | 40cm | Ø 1.6mm  | 5 Fr.  | <b>AC-019-002</b> | <b>AC-019-202</b>                                                                     | <b>AC-019-002-R</b> |
|  | Grasping forceps for stone street, double action | 40cm | Ø 1.2mm  | 4 Fr.  | <b>AC-007-001</b> | <b>AC-007-201</b>                                                                     | <b>AC-007-001-R</b> |
|                                                                                    |                                                  | 40cm | Ø 1.6mm  | 5 Fr.  | <b>AC-007-002</b> | <b>AC-007-202</b>                                                                     | <b>AC-007-002-R</b> |
|  | Grasping forceps for big stones, double action   | 40cm | Ø 1.2mm  | 4 Fr.  | <b>AC-009-001</b> | <b>AC-009-201</b>                                                                     | <b>AC-009-001-R</b> |
|                                                                                    |                                                  | 40cm | Ø 1.6mm  | 5 Fr.  | <b>AC-009-002</b> | <b>AC-009-202</b>                                                                     | <b>AC-009-002-R</b> |
|                                                                                    |                                                  | 40cm | Ø 2.3mm  | 7 Fr.  | <b>AC-009-004</b> | <b>AC-009-204</b>                                                                     | <b>AC-009-004-R</b> |
|  | 3-Prongs retriever                               | 40cm | Ø 1.0mm  | 3 Fr.  | <b>AC-012-001</b> |  |                     |
|                                                                                    |                                                  | 40cm | Ø 1.6mm  | 5 Fr.  | <b>AC-012-002</b> |                                                                                       |                     |
|                                                                                    |                                                  | 40cm | Ø 2.3mm  | 7 Fr.  | <b>AC-012-004</b> |                                                                                       |                     |

## 4.1 | FORCEPS, SCISSORS, INSTRUMENTS


### SCISSORS FLEXIBLE | SEMI-RIGID | RIGID



|                                                                                     | Description                               | WL   | Diameter | French | Flexible   | Semi-rigid | Rigid        |
|-------------------------------------------------------------------------------------|-------------------------------------------|------|----------|--------|------------|------------|--------------|
|    | Scissor, serrated flexible, double action | 40cm | Ø 2.3mm  | 7 Fr.  | AC-005-004 | AC-005-204 | AC-005-004-R |
|    | Scissor, single action                    | 40cm | Ø 1.0mm  | 3 Fr.  | AC-003-000 | AC-003-200 | AC-003-000-R |
|                                                                                     |                                           | 40cm | Ø 1.2mm  | 4 Fr.  | AC-003-001 | AC-003-201 | AC-003-001-R |
|                                                                                     |                                           | 40cm | Ø 1.6mm  | 5 Fr.  | AC-003-002 | AC-003-202 | AC-003-002-R |
|                                                                                     |                                           | 40cm | Ø 2.1mm  | 6 Fr.  | AC-003-003 | AC-003-203 | AC-003-003-R |
|                                                                                     |                                           | 40cm | Ø 2.3mm  | 7 Fr.  | AC-003-004 | AC-003-204 | AC-003-004-R |
|                                                                                     |                                           | 40cm | Ø 3.0mm  | 9 Fr.  | AC-003-005 | AC-003-205 | AC-003-005-R |
|   | Uretero slitting scissors, single action  | 40cm | Ø 1.6mm  | 5 Fr.  |            | AC-018-201 | AC-018-001-R |
|                                                                                     |                                           | 40cm | Ø 2.3mm  | 7 Fr.  |            | AC-018-202 | AC-018-002-R |
|   | Scissors pointed, single action           | 40cm | Ø 1.6mm  | 5 Fr.  |            | AC-003-302 |              |
|                                                                                     |                                           | 40cm | Ø 2.3mm  | 7 Fr.  |            | AC-003-304 |              |
|  | Hook scissors, single action              | 40cm | Ø 2.3mm  | 7 Fr.  | AC-004-004 | AC-004-204 | AC-004-004-R |
|                                                                                     |                                           | 40cm | Ø 3.0mm  | 9 Fr.  | AC-004-005 | AC-004-205 | AC-004-005-R |


## 4.1 | FORCEPS, SCISSORS, INSTRUMENTS

### INSTRUMENTS SEMI-RIGID



|                                                                                     | Description                                   | WL    | Diameter           | Order No.                              |
|-------------------------------------------------------------------------------------|-----------------------------------------------|-------|--------------------|----------------------------------------|
|    | Biopsy forceps oval, double action            | 340mm | Ø 1.6mm<br>Ø 2.3mm | <b>AC-001-502</b><br><b>AC-001-504</b> |
|                                                                                     |                                               | 410mm | Ø 1.6mm<br>Ø 2.3mm | <b>AC-001-512</b><br><b>AC-001-514</b> |
|   | Biopsy forceps oval, serrated, double action  | 340mm | Ø 1.6mm            | <b>AC-002-502</b>                      |
|                                                                                     |                                               | 410mm | Ø 1.6mm            | <b>AC-002-512</b>                      |
|  | Biopsy forceps oval, serrated, single action  | 340mm | Ø 1.6mm<br>Ø 2.3mm | <b>AC-010-502</b><br><b>AC-010-504</b> |
|                                                                                     |                                               | 410mm | Ø 1.6mm<br>Ø 2.3mm | <b>AC-010-512</b><br><b>AC-010-514</b> |
|  | Biopsy punch forceps, single action           | 340mm | Ø 1.6mm<br>Ø 2.3mm | <b>AC-014-502</b><br><b>AC-014-504</b> |
|                                                                                     |                                               | 410mm | Ø 1.6mm<br>Ø 2.3mm | <b>AC-014-512</b><br><b>AC-014-514</b> |
|  | Grasping forceps alligator jaw, double action | 340mm | Ø 1.6mm<br>Ø 2.3mm | <b>AC-006-502</b><br><b>AC-006-504</b> |
|                                                                                     |                                               | 410mm | Ø 1.6mm<br>Ø 2.3mm | <b>AC-006-512</b><br><b>AC-006-514</b> |
|  | Scissors, single action                       | 340mm | Ø 1.6mm<br>Ø 2.3mm | <b>AC-003-502</b><br><b>AC-003-504</b> |
|                                                                                     |                                               | 410mm | Ø 1.6mm<br>Ø 2.3mm | <b>AC-003-512</b><br><b>AC-003-514</b> |
|  | Scissors pointed, single action               | 340mm | Ø 1.6mm<br>Ø 2.3mm | <b>AC-003-602</b><br><b>AC-003-604</b> |
|                                                                                     |                                               | 410mm | Ø 1.6mm<br>Ø 2.3mm | <b>AC-003-612</b><br><b>AC-003-614</b> |

# NEW VERSION OF FORCEPS WITH LUER LOCK AND ROTARY WHEEL




## FEATURES

- 1 QUICK AND EASY LOCKING AND RELEASING OF THE HANDLE WITH CLICK LOCK
- 2 ROTARY WHEEL ENABLES PRECISE POSITIONING OF JAW SECTION
- 3 OVERLOAD PROTECTION AVAILABLE
- 4 SIDE-MOUNTED LUER LOCK

## 4.1 | FORCEPS, SCISSORS, INSTRUMENTS

### INSTRUMENTS SEMI-RIGID | DETACHABLE




| Ordering no.      | Description                              |
|-------------------|------------------------------------------|
| <b>AC-100-001</b> | <b>HANDLE FOR DETACHABLE INSTRUMENTS</b> |

|                                                                                    | Description                                      | Diameter | French | Complete            | Insert               |
|------------------------------------------------------------------------------------|--------------------------------------------------|----------|--------|---------------------|----------------------|
|   | Scissors, single action                          | Ø 1.0mm  | 3 Fr.  | <b>AC-003-200-T</b> | <b>AC-003-200-OT</b> |
|                                                                                    |                                                  | Ø 1.2mm  | 4 Fr.  | <b>AC-003-201-T</b> | <b>AC-003-201-OT</b> |
|                                                                                    |                                                  | Ø 1.6mm  | 5 Fr.  | <b>AC-003-202-T</b> | <b>AC-003-202-OT</b> |
|                                                                                    |                                                  | Ø 2.3mm  | 7 Fr.  | <b>AC-003-204-T</b> | <b>AC-003-204-OT</b> |
|                                                                                    |                                                  | Ø 3.0mm  | 9 Fr.  | <b>AC-003-205-T</b> | <b>AC-003-205-OT</b> |
|  | Scissors pointed, single action                  | Ø 1.6mm  | 5 Fr.  | <b>AC-003-612-T</b> | <b>AC-003-612-OT</b> |
|                                                                                    |                                                  | Ø 2.3mm  | 7 Fr.  | <b>AC-003-614-T</b> | <b>AC-003-614-OT</b> |
|  | Hook scissors, single action                     | Ø 2.3mm  | 7 Fr.  | <b>AC-004-204-T</b> | <b>AC-004-204-OT</b> |
|                                                                                    |                                                  | Ø 3.0mm  | 9 Fr.  | <b>AC-004-205-T</b> | <b>AC-004-205-OT</b> |
|  | Biopsy forceps oval, double action               | Ø 1.0mm  | 3 Fr.  | <b>AC-001-201-T</b> | <b>AC-001-201-OT</b> |
|                                                                                    |                                                  | Ø 1.2mm  | 4 Fr.  | <b>AC-001-206-T</b> | <b>AC-001-206-OT</b> |
|                                                                                    |                                                  | Ø 1.6mm  | 5 Fr.  | <b>AC-001-202-T</b> | <b>AC-001-202-OT</b> |
|                                                                                    |                                                  | Ø 2.3mm  | 7 Fr.  | <b>AC-001-203-T</b> | <b>AC-001-203-OT</b> |
|                                                                                    |                                                  | Ø 3.0mm  | 9 Fr.  | <b>AC-001-204-T</b> | <b>AC-001-204-OT</b> |
|  | Grasping forceps, alligator jaw, double action   | Ø 1.0mm  | 3 Fr.  | <b>AC-006-201-T</b> | <b>AC-006-201-OT</b> |
|                                                                                    |                                                  | Ø 1.2mm  | 4 Fr.  | <b>AC-006-203-T</b> | <b>AC-006-203-OT</b> |
|                                                                                    |                                                  | Ø 1.6mm  | 5 Fr.  | <b>AC-006-202-T</b> | <b>AC-006-202-OT</b> |
|                                                                                    |                                                  | Ø 2.3mm  | 7 Fr.  | <b>AC-006-204-T</b> | <b>AC-006-204-OT</b> |
|                                                                                    |                                                  | Ø 3.0mm  | 9 Fr.  | <b>AC-006-205-T</b> | <b>AC-006-205-OT</b> |
|  | Grasping forceps for stone street, double action | Ø 1.2mm  | 4 Fr.  | <b>AC-007-201-T</b> | <b>AC-007-201-OT</b> |
|                                                                                    |                                                  | Ø 1.6mm  | 5 Fr.  | <b>AC-007-202-T</b> | <b>AC-007-202-OT</b> |
|  | Grasping forceps for big stones, double action   | Ø 1.2mm  | 5 Fr.  | <b>AC-009-201-T</b> | <b>AC-009-201-OT</b> |
|                                                                                    |                                                  | Ø 1.6mm  | 5 Fr.  |                     | <b>AC-009-202-OT</b> |

## 4.1 | FORCEPS, SCISSORS, INSTRUMENTS

### INSTRUMENTS RIGID | DETACHABLE



| Ordering no.      | Description                              |
|-------------------|------------------------------------------|
| <b>AC-100-001</b> | <b>HANDLE FOR DETACHABLE INSTRUMENTS</b> |

|                                                                                     | Description                                      | Diameter | French | Complete            | Insert               |
|-------------------------------------------------------------------------------------|--------------------------------------------------|----------|--------|---------------------|----------------------|
|   | Scissors, single action                          | Ø 1.0mm  | 3 Fr.  | <b>AC-103-200-T</b> | <b>AC-103-200-OT</b> |
|                                                                                     |                                                  | Ø 1.2mm  | 4 Fr.  | <b>AC-103-201-T</b> | <b>AC-103-201-OT</b> |
|                                                                                     |                                                  | Ø 1.6mm  | 5 Fr.  | <b>AC-103-202-T</b> | <b>AC-103-202-OT</b> |
|                                                                                     |                                                  | Ø 2.3mm  | 7 Fr.  | <b>AC-103-204-T</b> | <b>AC-103-204-OT</b> |
|                                                                                     |                                                  | Ø 3.0mm  | 9 Fr.  | <b>AC-103-205-T</b> | <b>AC-103-205-OT</b> |
|  | Scissors pointed, single action                  | Ø 1.6mm  | 5 Fr.  | <b>AC-103-612-T</b> | <b>AC-103-612-OT</b> |
|                                                                                     |                                                  | Ø 2.3mm  | 7 Fr.  | <b>AC-103-614-T</b> | <b>AC-103-614-OT</b> |
|  | Hook scissors, single action                     | Ø 2.3mm  | 7 Fr.  | <b>AC-104-204-T</b> | <b>AC-104-204-OT</b> |
|                                                                                     |                                                  | Ø 3.0mm  | 9 Fr.  | <b>AC-104-205-T</b> | <b>AC-104-205-OT</b> |
|  | Biopsy forceps oval, double action               | Ø 1.0mm  | 3 Fr.  | <b>AC-101-201-T</b> | <b>AC-101-201-OT</b> |
|                                                                                     |                                                  | Ø 1.2mm  | 4 Fr.  | <b>AC-101-206-T</b> | <b>AC-101-206-OT</b> |
|                                                                                     |                                                  | Ø 1.6mm  | 5 Fr.  | <b>AC-101-202-T</b> | <b>AC-101-202-OT</b> |
|                                                                                     |                                                  | Ø 2.3mm  | 7 Fr.  | <b>AC-101-203-T</b> | <b>AC-101-203-OT</b> |
|                                                                                     |                                                  | Ø 3.0mm  | 9 Fr.  | <b>AC-101-204-T</b> | <b>AC-101-204-OT</b> |
|  | Grasping forceps, alligator jaw, double action   | Ø 1.0mm  | 3 Fr.  | <b>AC-106-201-T</b> | <b>AC-106-201-OT</b> |
|                                                                                     |                                                  | Ø 1.2mm  | 4 Fr.  | <b>AC-106-203-T</b> | <b>AC-106-203-OT</b> |
|                                                                                     |                                                  | Ø 1.6mm  | 5 Fr.  | <b>AC-106-202-T</b> | <b>AC-106-202-OT</b> |
|                                                                                     |                                                  | Ø 2.3mm  | 7 Fr.  | <b>AC-106-204-T</b> | <b>AC-106-204-OT</b> |
|                                                                                     |                                                  | Ø 3.0mm  | 9 Fr.  | <b>AC-106-205-T</b> | <b>AC-106-205-OT</b> |
|  | Grasping forceps for stone street, double action | Ø 1.2mm  | 4 Fr.  | <b>AC-107-201-T</b> | <b>AC-107-201-OT</b> |
|                                                                                     |                                                  | Ø 1.6mm  | 5 Fr.  | <b>AC-107-202-T</b> | <b>AC-107-202-OT</b> |
|  | Grasping forceps for big stones, double action   | Ø 1.2mm  | 5 Fr.  | <b>AC-109-201-T</b> | <b>AC-109-201-OT</b> |

## 4.2 | ELECTRODES & CONNECTING CABLES

### MONOPOLAR FLEXIBLE ELECTRODES



|                                                                                     | Description                          | French | Length | Ordering no. | Length | Ordering no. |
|-------------------------------------------------------------------------------------|--------------------------------------|--------|--------|--------------|--------|--------------|
|    | Loop electrode flexible, monopolar   | 5 Fr.  | 530mm  | AD-003-003   | 650mm  | AD-003-013   |
|                                                                                     |                                      | 7 Fr.  | 530mm  | AD-003-005   | 650mm  | AD-003-015   |
|                                                                                     |                                      | 9 Fr.  | 530mm  | AD-003-007   | 650mm  | AD-003-017   |
|                                                                                     |                                      | 10 Fr. | 530mm  | AD-003-008   | 650mm  | AD-003-018   |
|   | Needle electrode flexible, monopolar | 3 Fr.  | 530mm  | AD-002-001   | 650mm  | AD-002-011   |
|                                                                                     |                                      | 4 Fr.  | 530mm  | AD-002-002   | 650mm  | AD-002-012   |
|                                                                                     |                                      | 5 Fr.  | 530mm  | AD-002-003   | 650mm  | AD-002-013   |
|                                                                                     |                                      | 7 Fr.  | 530mm  | AD-002-005   | 650mm  | AD-002-015   |
|  | Button electrode flexible, monopolar | 3 Fr.  | 530mm  | AD-001-001   | 650mm  | AD-001-011   |
|                                                                                     |                                      | 4 Fr.  | 530mm  | AD-001-002   | 650mm  | AD-001-012   |
|                                                                                     |                                      | 5 Fr.  | 530mm  | AD-001-003   | 650mm  | AD-001-013   |
|                                                                                     |                                      | 6 Fr.  | 530mm  | AD-001-004   | 650mm  | AD-001-014   |
|                                                                                     |                                      | 7 Fr.  | 530mm  | AD-001-005   | 650mm  | AD-001-015   |
|                                                                                     |                                      | 8 Fr.  | 530mm  | AD-001-006   | 650mm  | AD-001-016   |
|                                                                                     |                                      | 9 Fr.  | 530mm  | AD-001-007   | 650mm  | AD-001-017   |
|                                                                                     |                                      | 10 Fr. | 530mm  | AD-001-008   | 650mm  | AD-001-018   |



## 4.2 | ELECTRODES & CONNECTING CABLES

### CONNECTING CABLES | MONOPOLAR



#### CONNECTION INSTRUMENT SIDE



| Unit Side                                                                          | Ordering no.      | Length | Diameter | Description                                                                  |
|------------------------------------------------------------------------------------|-------------------|--------|----------|------------------------------------------------------------------------------|
|   | <b>AD-010-001</b> | 3m     | Ø 4mm    | HF Monopolar connecting cable for <b>Berchtold, Martin, Aesculap</b>         |
|                                                                                    | <b>AD-010-006</b> | 5m     | Ø 4mm    |                                                                              |
|  | <b>AD-010-002</b> | 3m     | Ø 5mm    | HF Monopolar connecting cable for <b>Erbe, ACC/ICC, Storz</b>                |
|                                                                                    | <b>AD-010-007</b> | 5m     | Ø 5mm    |                                                                              |
|  | <b>AD-010-003</b> | 3m     | Ø 8mm    | HF Monopolar connecting cable for <b>STEMA HF-U Series, Valleylab/Conmed</b> |
|                                                                                    | <b>AD-010-008</b> | 5m     | Ø 8mm    |                                                                              |
|  | <b>AD-010-024</b> | 4.5m   | Ø 2mm    | HF Monopolar connecting cable <b>male</b> for <b>STEMA MASTER</b>            |

| Unit Side                                                                          | Ordering no.      | Length | Diameter | Description                                                                                           |
|------------------------------------------------------------------------------------|-------------------|--------|----------|-------------------------------------------------------------------------------------------------------|
|  | <b>AD-010-004</b> | 5m     | Ø 4mm    | HF Monopolar connecting cable <b>for units with 4mm banana plug</b><br><b>Not available in the EU</b> |
|                                                                                    | <b>AD-010-005</b> |        | Ø 4mm    |                                                                                                       |

## 4.2 | ELECTRODES & CONNECTING CABLES

### CONNECTING CABLES | MONOPOLAR



#### CONNECTION UNIT SIDE

| Instrument Side                                                                     | Ordering no.      | Length | Diameter             | Description                                                   |
|-------------------------------------------------------------------------------------|-------------------|--------|----------------------|---------------------------------------------------------------|
|  | <b>AD-010-025</b> | 4.5m   | 90° angled connector | HF Monopolar connecting cable for <b>STORZ resectoscope</b>   |
|  | <b>AD-010-026</b> | 4.5m   | Ø 3mm                | HF Monopolar connecting cable for <b>OLYMPUS resectoscope</b> |



#### RECOMMENDED PRODUCTS:

##### CI-001-030 STEMA HF-U PRO TOUCH

- For all monopolar & bipolar techniques
- with 7" touch screen display
- provided to perform the most advanced electrosurgery
- ADC / APC System

##### CI-001-040 DOUBLE PEDAL FOOTSWITCH

## 4.2 | ELECTRODES & CONNECTING CABLES

### BIPOLAR FLEXIBLE ELECTRODES



|                                                                                     | Description                                               | French | Length | Ordering no.      |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------|--------|--------|-------------------|
|    | Button electrode flexible, bipolar, flat plug             | 5 Fr.  | 360mm  | <b>AD-004-001</b> |
|    | Needle electrode angled, flexible, bipolar, flat plug     | 5 Fr.  | 360mm  | <b>AD-004-002</b> |
|   | Needle electrode 90° angled, flexible, bipolar, flat plug | 5 Fr.  | 360mm  | <b>AD-004-003</b> |
|  | Needle electrode straight, flexible, bipolar, flat plug   | 5 Fr.  | 360mm  | <b>AD-004-004</b> |

## 4.2 | ELECTRODES & CONNECTING CABLES

### CONNECTING CABLES | BIPOLAR



CONNECTION INSTRUMENT SIDE



| Connection Unit Side                                                                | Ordering no.      | Length | Description                                                                     |
|-------------------------------------------------------------------------------------|-------------------|--------|---------------------------------------------------------------------------------|
|  | <b>EK-010-010</b> | 3m     | Bipolar connecting cable, flat plug,<br>for <b>Berchtold, Martin, Aesculap</b>  |
|                                                                                     | <b>EK-010-011</b> | 5m     |                                                                                 |
|  | <b>EK-010-012</b> | 3m     | Bipolar connecting cable for, flat plug, <b>Erbe</b>                            |
|                                                                                     | <b>EK-010-013</b> | 5m     |                                                                                 |
|  | <b>EK-010-014</b> | 3m     | Bipolar connecting cable, flat plug,<br>for <b>STEMA, Valleylab/Conmed, EMC</b> |
|                                                                                     | <b>EK-010-015</b> | 5m     |                                                                                 |
|  | <b>EK-010-016</b> | 3m     | Bipolar connecting cable, flat plug, for units with<br><b>4mm banana plug</b>   |
|                                                                                     | <b>EK-010-017</b> | 5m     |                                                                                 |

## 4.3 | ACCESSORIES

### COHEN INTRA UTERINE PROBE




| Ordering no. | Description                                                      |
|--------------|------------------------------------------------------------------|
| AM-001-001   | COHEN intrauterine probe, <b>complete with 2 portio adapters</b> |
| AM-001-002   | COHEN intrauterine probe, <b>only</b>                            |
| AM-001-003   | Portio adapter, <b>small</b>                                     |
| AM-001-004   | Portio adapter, <b>large</b>                                     |



| Ordering no. | Description                                     |
|--------------|-------------------------------------------------|
| AM-004-001   | Ø 5mm myoma drill, working length <b>330mm</b>  |
| AM-004-002   | Ø 5mm myoma drill, working length <b>420mm</b>  |
| AM-004-004   | Ø 10mm myoma drill, working length <b>330mm</b> |
| AM-004-006   | Ø 10mm myoma drill, working length <b>420mm</b> |

## 4.3 | ACCESSORIES

### UTERUS MANIPULATOR



| Ordering no.      | Description                                         |
|-------------------|-----------------------------------------------------|
| <b>AM-002-001</b> | Uterus manipulator <b>complete</b> incl. 6 adaptors |
| <b>AM-002-002</b> | Uterus manipulator, <b>only</b>                     |

|                                                                                     | Ordering no.      | Description                                  |
|-------------------------------------------------------------------------------------|-------------------|----------------------------------------------|
|  | <b>AM-002-003</b> | Adaptor extra long, $\varnothing$ 10 x 100mm |
|  | <b>AM-002-004</b> | Adaptor long, $\varnothing$ 8 x 80mm         |
|  | <b>AM-002-005</b> | Adaptor thin, long, $\varnothing$ 3 x 55mm   |
|  | <b>AM-002-006</b> | Adaptor thin, short, $\varnothing$ 3 x 45mm  |
|  | <b>AM-002-007</b> | Adaptor with conus, <b>long</b>              |
|  | <b>AM-002-008</b> | Adaptor with conus, <b>short</b>             |

## 4.3 | ACCESSORIES

### UTERUS MANIPULATOR



AM-002-010



AM-002-011



AM-002-015



AM-002-016

| Ordering no. | Description                                                                                                           |
|--------------|-----------------------------------------------------------------------------------------------------------------------|
| AM-002-010   | Uterus manipulator set consisting of:<br><b>Basic instrument, 1 valve, 3 sealing foams and 11 adaptors (20-120mm)</b> |
| AM-002-011   | Adapter set for AM-002-010 consisting of 11 adaptors (20-120mm)                                                       |
| AM-002-012   | Ceramic cap small, Ø 24mm inner, Ø 32mm outer                                                                         |
| AM-002-013   | Ceramic cap medium, Ø 28mm inner, Ø 36mm outer                                                                        |
| AM-002-014   | Ceramic cap large, Ø 36mm inner, Ø 40mm outer                                                                         |
| AM-002-015   | Foam padded seal sterile, <b>3 pieces/box</b>                                                                         |
| AM-002-016   | Valve tube with hand lever, Ø 16mm                                                                                    |
| AM-002-018   | Basic instrument for uterus manipulator                                                                               |

## 4.3 | ACCESSORIES

### UTERUS MANIPULATOR | PALPATION PROBE



| Ordering no.      | Description                          |
|-------------------|--------------------------------------|
| <b>AM-002-020</b> | Uterus manipulator <b>complete</b>   |
| <b>AM-002-021</b> | Handle <b>rotatable</b>              |
| <b>AM-002-022</b> | Sheath with fixation nut             |
| <b>AM-002-023</b> | Inner sheath                         |
| <b>AM-002-024</b> | Sliding and fixation sheath          |
| <b>AM-002-025</b> | Tube <b>35mm</b>                     |
| <b>AM-002-026</b> | Tube <b>40mm</b>                     |
| <b>AM-002-027</b> | Tube <b>45mm</b>                     |
| <b>AM-002-028</b> | Cervical tenaculum forceps           |
| <b>AM-002-029</b> | Intrauterin adapter, <b>4 x 40mm</b> |
| <b>AM-002-030</b> | Intrauterin adapter, <b>5 x 50mm</b> |
| <b>AM-002-031</b> | Intrauterin adapter, <b>6 x 60mm</b> |
| <b>AM-002-032</b> | Intrauterin adapter, <b>7 x 70mm</b> |




| Ordering no.      | Description                                                     |
|-------------------|-----------------------------------------------------------------|
| <b>AN-002-020</b> | Palpation probe with scaling, semi-rigid, <b>5 Fr., WL 34cm</b> |



## 4.3 | ACCESSORIES

### TENACULUM FORCEPS



AJ-005-001



AJ-005-002

| Ordering no. | Description                                             |
|--------------|---------------------------------------------------------|
| AJ-005-001   | Barret, tenaculum forceps, <b>18cm</b>                  |
| AJ-005-002   | Pozzi, tenaculum forceps, strong pattern, <b>25.5cm</b> |

## 4.3 | ACCESSORIES

### DILATATOR BOUGIES



VAN BUREN

DITTEL,  
STRAIGHTDITTEL,  
CURVED

| French | Ordering no.<br>VAN BUREN | Ordering no.<br>DITTEL straight | Ordering no.<br>DITTEL curved |
|--------|---------------------------|---------------------------------|-------------------------------|
| 8 Fr.  | AJ-003-001                | AJ-003-021                      | AJ-003-050                    |
| 9 Fr.  |                           | AJ-003-022                      | AJ-003-051                    |
| 10 Fr. | AJ-003-002                | AJ-003-023                      | AJ-003-052                    |
| 11 Fr. |                           | AJ-003-024                      | AJ-003-053                    |
| 12 Fr. | AJ-003-003                | AJ-003-025                      | AJ-003-054                    |
| 13 Fr. |                           | AJ-003-026                      | AJ-003-055                    |
| 14 Fr. | AJ-003-004                | AJ-003-027                      | AJ-003-056                    |
| 15 Fr. |                           | AJ-003-028                      | AJ-003-057                    |
| 16 Fr. | AJ-003-005                | AJ-003-029                      | AJ-003-058                    |
| 17 Fr. |                           | AJ-003-030                      | AJ-003-059                    |
| 18 Fr. | AJ-003-006                | AJ-003-031                      | AJ-003-060                    |
| 19 Fr. |                           | AJ-003-032                      | AJ-003-061                    |
| 20 Fr. | AJ-003-007                | AJ-003-033                      | AJ-003-062                    |
| 21 Fr. |                           | AJ-003-034                      | AJ-003-063                    |
| 22 Fr. | AJ-003-008                | AJ-003-035                      | AJ-003-064                    |
| 23 Fr. |                           | AJ-003-036                      | AJ-003-065                    |
| 24 Fr. | AJ-003-009                | AJ-003-037                      | AJ-003-066                    |
| 25 Fr. |                           | AJ-003-038                      | AJ-003-067                    |
| 26 Fr. | AJ-003-010                | AJ-003-039                      | AJ-003-068                    |
| 27 Fr. |                           | AJ-003-040                      | AJ-003-069                    |
| 28 Fr. | AJ-003-011                | AJ-003-041                      | AJ-003-070                    |
| 29 Fr. |                           | AJ-003-042                      | AJ-003-071                    |
| 30 Fr. | AJ-003-012                | AJ-003-043                      | AJ-003-072                    |
| 31 Fr. |                           |                                 |                               |
| 32 Fr. | AJ-003-013                | AJ-003-044                      | AJ-003-074                    |
| 33 Fr. |                           |                                 |                               |
| 34 Fr. | AJ-003-014                | AJ-003-045                      |                               |
| 35 Fr. |                           |                                 |                               |
| 36 Fr. | AJ-003-015                | AJ-003-046                      |                               |
| 37 Fr. |                           |                                 |                               |
| 38 Fr. | AJ-003-016                |                                 |                               |

## 4.3 | ACCESSORIES

### DILATATOR BOUGIES



GUYON


| French |            | Ordering no. GUYON |            |
|--------|------------|--------------------|------------|
| 8 Fr.  | AJ-003-170 | 25 Fr.             | AJ-003-187 |
| 9 Fr.  | AJ-003-171 | 26 Fr.             | AJ-003-188 |
| 10 Fr. | AJ-003-172 | 27 Fr.             | AJ-003-189 |
| 11 Fr. | AJ-003-173 | 28 Fr.             | AJ-003-190 |
| 12 Fr. | AJ-003-174 | 29 Fr.             | AJ-003-191 |
| 13 Fr. | AJ-003-175 | 30 Fr.             | AJ-003-192 |
| 14 Fr. | AJ-003-176 | 31 Fr.             | AJ-003-193 |
| 15 Fr. | AJ-003-177 | 32 Fr.             | AJ-003-194 |
| 16 Fr. | AJ-003-178 | 33 Fr.             | AJ-003-195 |
| 17 Fr. | AJ-003-179 | 34 Fr.             | AJ-003-196 |
| 18 Fr. | AJ-003-180 | 35 Fr.             | AJ-003-197 |
| 19 Fr. | AJ-003-181 | 36 Fr.             | AJ-003-198 |
| 20 Fr. | AJ-003-182 | 37 Fr.             | AJ-003-199 |
| 21 Fr. | AJ-003-183 | 38 Fr.             | AJ-003-200 |
| 22 Fr. | AJ-003-184 | 39 Fr.             | AJ-003-201 |
| 23 Fr. | AJ-003-185 | 40 Fr.             | AJ-003-202 |
| 24 Fr. | AJ-003-186 | 50 Fr.             | AJ-003-203 |

HEGAR

| French    |            | Ordering no. HEGAR |  |
|-----------|------------|--------------------|--|
| 3/4 Fr.   | AJ-003-151 |                    |  |
| 5/6 Fr.   | AJ-003-152 |                    |  |
| 7/8 Fr.   | AJ-003-153 |                    |  |
| 9/10 Fr.  | AJ-003-154 |                    |  |
| 11/12 Fr. | AJ-003-155 |                    |  |
| 13/14 Fr. | AJ-003-156 |                    |  |
| 15/16 Fr. | AJ-003-157 |                    |  |
| 17/18 Fr. | AJ-003-158 |                    |  |
|           |            |                    |  |
|           |            |                    |  |
|           |            |                    |  |
|           |            |                    |  |
|           |            |                    |  |
|           |            |                    |  |
|           |            |                    |  |
|           |            |                    |  |

## 4.3 | ACCESSORIES

### METAL CATHETERS



FOR MEN


FOR WOMEN

COXETER  
PROSTATA

| French | Ordering no. | Ordering no. | Ordering no. |
|--------|--------------|--------------|--------------|
| 6 Fr.  |              | AJ-003-100   |              |
| 7 Fr.  | AJ-003-081   |              |              |
| 8 Fr.  | AJ-003-082   | AJ-003-102   | AJ-003-120   |
| 9 Fr.  |              |              |              |
| 10 Fr. | AJ-003-084   | AJ-003-104   | AJ-003-121   |
| 11 Fr. |              |              |              |
| 12 Fr. | AJ-003-086   | AJ-003-106   | AJ-003-122   |
| 13 Fr. |              |              |              |
| 14 Fr. | AJ-003-088   | AJ-003-108   | AJ-003-123   |
| 15 Fr. |              |              |              |
| 16 Fr. | AJ-003-090   | AJ-003-110   | AJ-003-124   |
| 17 Fr. |              |              |              |
| 18 Fr. | AJ-003-092   | AJ-003-112   | AJ-003-125   |
| 19 Fr. |              |              |              |
| 20 Fr. | AJ-003-094   | AJ-003-114   | AJ-003-126   |
| 21 Fr. |              |              |              |
| 22 Fr. | AJ-003-096   | AJ-003-116   | AJ-003-127   |
| 23 Fr. |              |              |              |
| 24 Fr. | AJ-003-098   | AJ-003-118   | AJ-003-128   |

## 4.3 | ACCESSORIES

### UTERINE PROBES



| Ordering no. | Description                                  |
|--------------|----------------------------------------------|
| AJ-006-001   | Sims uterine probe <b>rigid</b> , 32cm       |
| AJ-006-002   | Sims uterine probe <b>malleable</b> , 32cm   |
| AJ-006-003   | Martin uterine probe <b>rigid</b> , 32cm     |
| AJ-006-004   | Martin uterine probe <b>malleable</b> , 32cm |

### BIOPSY CURETTES



| Ordering no. | Description                                   |
|--------------|-----------------------------------------------|
| AJ-006-005   | Novak endometrium biopsy curette, Ø 1mm, 23cm |
| AJ-006-006   | Novak endometrium biopsy curette, Ø 2mm, 23cm |
| AJ-006-007   | Novak endometrium biopsy curette, Ø 3mm, 23cm |
| AJ-006-008   | Novak endometrium biopsy curette, Ø 4mm, 23cm |

## 4.3 | ACCESSORIES

### HEGAR UTERINE DILATATOR




| Ordering no.      | Description                                                                        |
|-------------------|------------------------------------------------------------------------------------|
| <b>AJ-003-078</b> | Hegar uterine dilatation set ( <b>8 pieces</b> ), double ended, Ø 3/4mm - 7/18mm   |
| <b>AJ-003-079</b> | Hegar uterine dilatation set ( <b>13 pieces</b> ), double ended, Ø 1/2mm - 25/26mm |

### SECRETION SPOON




| Ordering no.      | Description                     |
|-------------------|---------------------------------|
| <b>AJ-006-012</b> | Guttman secretion spoon, 26.5cm |

### RETRACTOR




| Ordering no.      | Description                                                          |
|-------------------|----------------------------------------------------------------------|
| <b>AM-002-040</b> | LOWSLEY retractor, <b>straight</b> , Ø 21.5Fr., working length 280mm |
| <b>AM-002-041</b> | LOWSLEY retractor, <b>curved</b> , Ø 21.5Fr., working length 280mm   |

## 4.3 | ACCESSORIES | STERILIZATION



| Ordering no.      | Dimension<br>(l x w x h) | Description                                  |
|-------------------|--------------------------|----------------------------------------------|
| <b>FA-010-001</b> | <b>290 x 80 x 52mm</b>   | Wire basket for <b>1 endoscope</b> with lid  |
| <b>FA-010-002</b> | <b>460 x 80 x 52mm</b>   |                                              |
| <b>FA-010-004</b> | <b>290 x 80 x 52mm</b>   | Wire basket for <b>2 endoscopes</b> with lid |
| <b>FA-010-005</b> | <b>460 x 80 x 52mm</b>   |                                              |



| Ordering no.      | Dimension<br>(l x w x h) | Description                                                               |
|-------------------|--------------------------|---------------------------------------------------------------------------|
| <b>FA-010-006</b> | <b>670 x 80 x 52mm</b>   | Wire basket for <b>1 endoscope</b> and <b>1 cold light cable</b> with lid |

## 4.3 | STERILIZATION



| Ordering no. | Colour | Dimension<br>(l x w x h) | Description                                                        |
|--------------|--------|--------------------------|--------------------------------------------------------------------|
| FA-001-001   | Silver | 450 x 70 x 70mm          | Endoscope container, lid and bottom perforated<br>(without filter) |
| FA-001-001-R | Red    |                          |                                                                    |
| FA-001-001-G | Green  |                          |                                                                    |
| FA-002-001   | Silver | 250 x 60 x 50mm          | Endoscope container, lid and bottom perforated<br>(without filter) |
| FA-002-001-R | Red    |                          |                                                                    |
| FA-002-001-G | Green  |                          |                                                                    |



| Ordering no. | Description                                                  |
|--------------|--------------------------------------------------------------|
| FA-001-000   | Protection tube for endoscopes, plastic, 10 pcs, <b>19cm</b> |



## 5. SALINE RESECTOSCOPY

5.1 | 24/26 FRENCH

5.2 | 19/22 FRENCH (SLIM)

5.3 | 17.5 FRENCH (MINI)

5.4 | CONNECTING CABLES

# 5. SALINE RESECTOSCOPY

24/26 FR.  
28.5 FR.

19/22 FR.

17.5 FR.


## COMPLETE RESECTOSCOPE




### 1 CYSTOSCOPES



### 2 WORKING ELEMENTS



### 3 SHEATHS & OBTURATORS



### 4 ELECTRODES



## 5.1 | SALINE RESECTOSCOPY

### CYSTOSCOPES AUTOCLAVABLE Ø 4.0MM




| Viewing angle                                                                            | Endoscope                  | Diameter | Working length | Ordering no. |
|------------------------------------------------------------------------------------------|----------------------------|----------|----------------|--------------|
| <br>12° | HD Cystoscope autoclavable | Ø 4.0mm  | 302mm          | AA-002-012   |

| Viewing angle                                                                              | Endoscope                        | Diameter | Working length | Ordering no. |
|--------------------------------------------------------------------------------------------|----------------------------------|----------|----------------|--------------|
| <br>30° | Cystoscope autoclavable          | Ø 4.0mm  | 302mm          | AA-002-003   |
|                                                                                            | HD Cystoscope autoclavable       | Ø 4.0mm  | 302mm          | AA-002-013   |
|                                                                                            | ULTRA HD Cystoscope autoclavable | Ø 4.0mm  | 302mm          | AA-002-023   |

## 5.1 | SALINE RESECTOSCOPY

### CONTINUOUS FLOW SHEATHS | ROTATABLE



|                                                                                     | Ordering no.      | Description                                                             | French          | Colour code |
|-------------------------------------------------------------------------------------|-------------------|-------------------------------------------------------------------------|-----------------|-------------|
|    | <b>AH-001-302</b> | Outer sheath with perforation, <b>rotatable</b>                         | <b>26 Fr.</b>   | Yellow      |
|    | <b>AH-001-312</b> | Outer sheath with perforation, <b>rotatable</b>                         | <b>28.5 Fr.</b> | Black       |
|    | <b>AH-001-304</b> | Outer sheath with upper/lower perforation, <b>rotatable</b>             | <b>26 Fr.</b>   | Yellow      |
|   | <b>AH-001-306</b> | Outer sheath with surrounding fenestrated perforation, <b>rotatable</b> | <b>26 Fr.</b>   | Yellow      |
|  | <b>AH-001-303</b> | Outer sheath with upper fenestrated perforation, <b>rotatable</b>       | <b>26 Fr.</b>   | Yellow      |
|  | <b>AH-001-305</b> | Outer sheath with upper fenestrated perforation, <b>rotatable</b>       | <b>26 Fr.</b>   | Yellow      |



| Ordering no.      | Description                                           | French          | Colour code |
|-------------------|-------------------------------------------------------|-----------------|-------------|
| <b>AH-001-301</b> | Inner sheath with ceramic isolation, <b>rotatable</b> | <b>26 Fr.</b>   | Yellow      |
| <b>AH-001-311</b> | Inner sheath with ceramic isolation, <b>rotatable</b> | <b>28.5 Fr.</b> | Black       |

## 5.1 | SALINE RESECTOSCOPY OBTURATORS



| Ordering no. | Description                                        | French | For outer sheath | Colour code |
|--------------|----------------------------------------------------|--------|------------------|-------------|
| AH-001-011   | <b>Standard</b> Obturator for resectoscope sheaths | 24 Fr. | 26 Fr.           | Yellow      |
| AH-001-012   | <b>Standard</b> Obturator for resectoscope sheaths | 27 Fr. | 28.5 Fr.         | Black       |



| Ordering no. | Description                                          | French | For outer sheath | Colour code |
|--------------|------------------------------------------------------|--------|------------------|-------------|
| AH-001-015   | <b>Deflecting</b> Obturator for resectoscope sheaths | 24 Fr. | 26 Fr.           | Yellow      |
| AH-001-016   | <b>Deflecting</b> Obturator for resectoscope sheaths | 27 Fr. | 28.5 Fr.         | Black       |



| Ordering no. | Description                                       | French | For outer sheath | Colour code |
|--------------|---------------------------------------------------|--------|------------------|-------------|
| AH-001-013   | <b>Viewing</b> Obturator for resectoscope sheaths | 24 Fr. | 26 Fr.           | Yellow      |



| Ordering no. | Description                                                                            | French | For outer sheath | Colour code |
|--------------|----------------------------------------------------------------------------------------|--------|------------------|-------------|
| AH-001-017   | <b>Viewing</b> Obturator for resectoscope sheaths with 5 Fr. <b>Instrument channel</b> | 24 Fr. | 26 Fr.           | Yellow      |
| AH-001-018   | <b>Viewing</b> Obturator for resectoscope sheaths with 5 Fr. <b>Instrument channel</b> | 27 Fr. | 28.5 Fr.         | Black       |

# 5.1 | SALINE RESECTOSCOPY WORKING ELEMENTS



| Ordering no. | Description                                                                                                                      |
|--------------|----------------------------------------------------------------------------------------------------------------------------------|
| AG-001-011   | SALINE working element used for resectoscopes, <b>24 Fr.</b> , <b>passive cutting</b> , handle closed, cutting spring assistance |



| Ordering no. | Description                                                                                                                     |
|--------------|---------------------------------------------------------------------------------------------------------------------------------|
| AG-001-012   | SALINE working element used for resectoscopes, <b>24 Fr.</b> , <b>active cutting</b> , handle closed, cutting spring assistance |

## 5.1 | SALINE RESECTOSCOPY

### ELECTRODES FOR SALINE RESECTOSCOPY



|                                                                                     | Ordering no.      | French        | Description                                             |
|-------------------------------------------------------------------------------------|-------------------|---------------|---------------------------------------------------------|
|    | <b>AE-001-030</b> | <b>24 Fr.</b> | Cutting loop electrode <b>30° angled, strong cut</b>    |
|    | <b>AE-001-031</b> | <b>24 Fr.</b> | Cutting loop electrode <b>30° angled</b>                |
|    | <b>AE-001-033</b> | <b>24 Fr.</b> | Cutting loop electrode <b>90 ° angled</b>               |
|  | <b>AE-001-048</b> | <b>24 Fr.</b> | Cutting loop electrode <b>90 ° angled, for 0° Scope</b> |
|  | <b>AE-001-036</b> | <b>24 Fr.</b> | Cutting loop electrode <b>straight</b>                  |
|  | <b>AE-001-032</b> | <b>24 Fr.</b> | Knife electrode                                         |
|  | <b>AE-001-041</b> | <b>24 Fr.</b> | Band electrode, <b>30° angled</b>                       |
|  | <b>AE-001-038</b> | <b>24 Fr.</b> | Vaporisation spike rollerball electrode <b>Ø 3mm</b>    |
|  | <b>AE-001-037</b> | <b>24 Fr.</b> | Vaporisation rollerball electrode <b>Ø 3mm</b>          |
|  | <b>AE-001-043</b> | <b>24 Fr.</b> | Vaporisation rollerball electrode <b>Ø 5mm</b>          |
|  | <b>AE-001-040</b> | <b>24 Fr.</b> | Roller electrode smooth <b>Ø 3mm</b>                    |
|  | <b>AE-001-042</b> | <b>24 Fr.</b> | Roller electrode smooth <b>Ø 5mm</b>                    |


## 5.1 | SALINE RESECTOSCOPY

### ELECTRODES FOR SALINE RESECTOSCOPY



|                                                                                     | Ordering no.      | French        | Description                                  |
|-------------------------------------------------------------------------------------|-------------------|---------------|----------------------------------------------|
|    | <b>AE-001-039</b> | <b>24 Fr.</b> | Conical electrode                            |
|    | <b>AE-001-044</b> | <b>24 Fr.</b> | Conical electrode „Mushroom“                 |
|    | <b>AE-001-045</b> | <b>24 Fr.</b> | Conical electrode „Mushroom“, <b>convex</b>  |
|  | <b>AE-001-046</b> | <b>24 Fr.</b> | Conical electrode „Mushroom“, <b>concav</b>  |
|  | <b>AE-001-047</b> | <b>24 Fr.</b> | Conical electrode „Mushroom“, <b>grooved</b> |

|                                                                                     | Ordering no.      | French        | Description                               |
|-------------------------------------------------------------------------------------|-------------------|---------------|-------------------------------------------|
|  | <b>AE-002-031</b> | <b>27 Fr.</b> | Cutting loop electrode, <b>30° angled</b> |



| Ordering no.      | Description                                                                                 |
|-------------------|---------------------------------------------------------------------------------------------|
| <b>AJ-010-006</b> | Protection sleeve made of steel for sterilization and storage of electrodes and cold knives |



## 5.2 | SALINE SLIM RESECTOSCOPY

### CYSTOSCOPES AUTOCLAVABLE Ø 2.9MM



| Viewing angle                                                                            | Endoscope                  | Diameter | Working length | Ordering no. |
|------------------------------------------------------------------------------------------|----------------------------|----------|----------------|--------------|
| <br>12° | Cystoscope autoclavable    | Ø 2.9mm  | 302mm          | AA-003-004   |
|                                                                                          | HD Cystoscope autoclavable | Ø 2.9mm  | 302mm          | AA-003-014   |

| Viewing angle                                                                              | Endoscope                  | Diameter | Working length | Ordering no. |
|--------------------------------------------------------------------------------------------|----------------------------|----------|----------------|--------------|
| <br>30° | Cystoscope autoclavable    | Ø 2.9mm  | 302mm          | AA-003-002   |
|                                                                                            | HD Cystoscope autoclavable | Ø 2.9mm  | 302mm          | AA-003-012   |

## 5.2 | SALINE SLIM RESECTOSCOPY SHEATHS



| Ordering no. | Description                                            |
|--------------|--------------------------------------------------------|
| AH-002-302   | Outer sheath 22 Fr. <b>with perforation, rotatable</b> |



| Ordering no. | Description                                                  |
|--------------|--------------------------------------------------------------|
| AH-002-301   | Inner sheath 22 Fr. <b>with ceramic isolation, rotatable</b> |

## 5.2 | SALINE SLIM RESECTOSCOPY OBTURATORS



| Ordering no. | Description                                                      |
|--------------|------------------------------------------------------------------|
| AH-002-011   | <b>Standard obturator</b> for resectoscope sheaths <b>19 Fr.</b> |



| Ordering no. | Description                                                     |
|--------------|-----------------------------------------------------------------|
| AH-002-013   | <b>Viewing Obturator</b> for resectoscope sheaths <b>19 Fr.</b> |

## 5.2 | SALINE SLIM RESECTOSCOPY

### WORKING ELEMENTS



| Ordering no. | Description                                                                                                              |
|--------------|--------------------------------------------------------------------------------------------------------------------------|
| AG-002-011   | SALINE working element used for resectoscopes, <b>19 Fr., passive cutting</b> , handle closed, cutting spring assistance |



| Ordering no. | Description                                                                                                             |
|--------------|-------------------------------------------------------------------------------------------------------------------------|
| AG-002-012   | SALINE working element used for resectoscopes, <b>19 Fr., active cutting</b> , handle closed, cutting spring assistance |

## 5.2 | SALINE SLIM RESECTOSCOPY

### ELECTRODES FOR SALINE RESECTOSCOPY



|                                                                                     | Ordering no.      | French        | Description                              |
|-------------------------------------------------------------------------------------|-------------------|---------------|------------------------------------------|
|    | <b>AE-003-031</b> | <b>19 Fr.</b> | Cutting loop electrode <b>30° angled</b> |
|    | <b>AE-003-032</b> | <b>19 Fr.</b> | Cutting loop electrode <b>90° angled</b> |
|    | <b>AE-003-033</b> | <b>19 Fr.</b> | Cutting loop electrode <b>straight</b>   |
|  | <b>AE-003-034</b> | <b>19 Fr.</b> | Rollerball electrode <b>Ø 3mm</b>        |
|  | <b>AE-003-040</b> | <b>19 Fr.</b> | Roller electrode <b>Ø 3mm</b>            |
|  | <b>AE-003-036</b> | <b>19 Fr.</b> | Knife electrode                          |



| Ordering no.      | Description                                                                                 |
|-------------------|---------------------------------------------------------------------------------------------|
| <b>AJ-010-006</b> | Protection sleeve made of steel for sterilization and storage of electrodes and cold knives |



17.5 FR.


5.3 | SALINE MINI RESECTOSCOPY  
CYSTOSCOPES



| Viewing angle                                                                            | Endoscope                  | Diameter | Working length | Ordering no. |
|------------------------------------------------------------------------------------------|----------------------------|----------|----------------|--------------|
| <br>30° | Cystoscope autoclavable    | Ø 2.9mm  | 302mm          | AA-003-002   |
|                                                                                          | HD Cystoscope autoclavable | Ø 2.9mm  | 302mm          | AA-003-012   |

## 5.3 | SALINE MINI RESECTOSCOPY

### WORKING ELEMENTS



| Ordering no.      | Description                                                                                                                        |
|-------------------|------------------------------------------------------------------------------------------------------------------------------------|
| <b>AG-002-015</b> | SALINE working element used for resectoscopes, <b>18.5 Fr.</b> , <b>passive cutting</b> , handle closed, cutting spring assistance |



| Ordering no.      | Description                                                                                                                       |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| <b>AG-002-013</b> | SALINE working element used for resectoscopes, <b>18.5 Fr.</b> , <b>active cutting</b> , handle closed, cutting spring assistance |

## 5.3 | SALINE MINI RESECTOSCOPY SHEATHS



| Ordering no.      | Description                                   |
|-------------------|-----------------------------------------------|
| <b>AH-002-213</b> | Outer sheath 18.5 Fr. <b>with perforation</b> |



| Ordering no.      | Description                                                    |
|-------------------|----------------------------------------------------------------|
| <b>AH-002-211</b> | Inner sheath 17.5 Fr. <b>with ceramic isolation, rotatable</b> |



| Ordering no.      | Description                                                                                                      |
|-------------------|------------------------------------------------------------------------------------------------------------------|
| <b>AH-002-216</b> | Inner sheath 17.5 Fr. <b>without ceramic isolation</b> , with laser channel up to 5 Fr., for use with AH-001-213 |



### 5.3 | SALINE MINI RESECTOSCOPY OBTURATORS



| Ordering no. | Description                 |
|--------------|-----------------------------|
| AH-002-015   | Standard obturator 17.5 Fr. |



| Ordering no. | Description                                              |
|--------------|----------------------------------------------------------|
| AH-002-016   | Viewing Obturator 17.5 Fr. with 5 fr. instrument channel |

## 5.3 | SALINE MINI RESECTOSCOPY

### ELECTRODES FOR BIPOLAR RESECTOSCOPY




|                                                                                     | Ordering no.      | French          | Description                              |
|-------------------------------------------------------------------------------------|-------------------|-----------------|------------------------------------------|
|    | <b>AE-003-041</b> | <b>17.5 Fr.</b> | Cutting loop electrode <b>30° angled</b> |
|    | <b>AE-003-043</b> | <b>17.5 Fr.</b> | Cutting loop electrode <b>straight</b>   |
|    | <b>AE-003-044</b> | <b>17.5 Fr.</b> | Rollerball electrode <b>Ø 2mm</b>        |
|  | <b>AE-003-045</b> | <b>17.5 Fr.</b> | Roller electrode <b>Ø 2mm</b>            |
|  | <b>AE-003-046</b> | <b>17.5 Fr.</b> | Knife electrode                          |



| Ordering no.      | Description                                                                                 |
|-------------------|---------------------------------------------------------------------------------------------|
| <b>AJ-010-006</b> | Protection sleeve made of steel for sterilization and storage of electrodes and cold knives |

## 5.4 | SALINE RESECTOSCOPY

### CONNECTING CABLES | BIPOLAR



#### CONNECTION INSTRUMENT SIDE



| Connection Unit Side                                                                | Ordering no.      | Length | Description                                                              |
|-------------------------------------------------------------------------------------|-------------------|--------|--------------------------------------------------------------------------|
|  | <b>AD-010-011</b> | 4.5m   | Bipolar resectoscope connecting cable for <b>Martin Maxium</b>           |
|  | <b>AD-010-012</b> | 4.5m   | Bipolar resectoscope connecting cable for <b>Erbe Vio</b>                |
|  | <b>AD-010-013</b> | 4.5m   | Bipolar resectoscope connecting cable for <b>Valleylab, Lamidey, EMC</b> |
|  | <b>AD-010-015</b> | 4.5m   | Bipolar resectoscope connecting cable for <b>Erbe Vio 3</b>              |

## 6. HYSTEROSCOPY SETS

6.1 | DIAGNOSTIC HYSTEROSCOPY SET

6.2 | OPERATING HYSTEROSCOPY SET

6.3 | OUT-PATIENT ATRAUMATIC OP-HYSTEROSCOPY SET

6.4 | SALINE RESECTOSCOPY SET

# 6.1 CONTINUOUS-FLOW DIAGNOSTIC HYSTEROSCOPY SET

## 1 HD CYSTOSCOPE

Ø 4.0mm, 30° viewing angle, autoclavable, 302mm working length

1 AA-002-013



## 2 HYSTEROSCOPY INNER SHEATH

Inner-Diagnostic-Continuous-Flow sheath for 4mm / 30° scope, Ø 5.4mm

2 AN-003-102



## 3 HYSTEROSCOPY OUTER SHEATH

Outer-Diagnostic-Continuous-Flow sheath for 4mm scopes, Ø 6.5mm

3 AN-003-103



## 4 OBTURATOR

Used for Diagnostic-Continuous-Flow hysteroscope sheaths for 4mm scopes

4 AN-003-104



## 6.2 CONTINUOUS-FLOW OPERATING HYSTEROSCOPY SET

### 1 HD CYSTOSCOPE

Ø 4.0mm, 30° viewing angle, autoclavable, 302mm working length

1 AA-002-013



### 2 HYSTEROSCOPY OUTER SHEATH

Outer-Diagnostic-Continuous-Flow sheath for 4mm / 30° scope, Ø 6.5mm

2 AN-003-302



### 3 HYSTEROSCOPY INNER SHEATH

Inner-Diagnostic-Continuous-Flow sheath for 4mm scopes, Ø 5.4mm

3 AH-003-301



### 4 BIOPSY FORCEPS

oval, flexible, double action, Ø 2.3mm, WL 40cm

4 AC-001-004



### 5 SCISSOR FLEXIBLE

single action, Ø 2.3mm, WL 40cm

5 AC-003-004



### 6 GRASPING FORCEPS

alligator jaw flexible, double action, Ø 2.3mm, WL 40cm

6 AC-006-004



### 7 BUTTON ELECTRODE FLEXIBLE

7Fr., monopolar, length 530mm

7 AD-001-005



### 8 NEEDLE ELECTRODE FLEXIBLE

7Fr., monopolar, length 530mm

8 AD-002-005



## 6.3 OUT-PATIENT ATRAUMATIC OPERATING HYSTEROSCOPY SET

### 1 HD HYSTEROSCOPE

Ø 2.9mm, 30° viewing angle, autoclavable, 302mm working length

1 AA-003-012



### 2 HYSTEROSCOPY OUTER SHEATH

Outer-Diagnostic-Continuous-Flow sheath for 2.9mm / 30° scope, Ø 6.0mm

2 AN-002-302



### 3 HYSTEROSCOPY INNER SHEATH

Inner-Diagnostic-Continuous-Flow sheath for 2.9mm scopes, Ø 5.3mm

3 AH-002-301



### 4 BIOPSY FORCEPS

oval, flexible, double action, Ø 1.6mm, WL 40cm

4 AC-001-002



### 5 SCISSOR FLEXIBLE

single action, Ø 1.6mm, WL 40cm

5 AC-003-002



### 6 GRASPING FORCEPS

alligator jaw flexible, double action, Ø 1.6mm, WL 40cm

6 AC-006-002



### 7 BUTTON ELECTRODE FLEXIBLE

5Fr., monopolar, length 530mm

7 AD-001-003



### 8 NEEDLE ELECTRODE FLEXIBLE

5Fr., monopolar, length 530mm

8 AD-002-003



## 6.4 SALINE MINI RESECTOSCOPY SET

### 1 HD CYSTOSCOPE

Ø 2.9mm, 30° viewing angle, autoclavable,  
WL 302mm

1 AA-003-012



### 2 WORKING ELEMENT

for resectoscopes, 18.5 FR., passive cutting,  
handle closed, cutting spring assistance

2 AG-002-015



### 3 OUTER SHEATH

18.5 Fr. with perforation

3 AH-002-213



### 4 INNER SHEATH

17.5 Fr. with ceramic isolation, rotatable

4 AH-002-211



### 5 STANDARD OBTURATOR

17.5 Fr.

5 AH-002-015



### 6 CUTTING LOOP ELECTRODE

30° angled, 17.5 Fr.

6 AE-003-041



### 7 ROLLERBALL ELECTRODE 3 MM

17.5 Fr.

7 AE-003-044



### 8 KNIFE ELECTRODE

17.5 Fr.

8 AE-003-046





ARE YOU INTERESTED  
IN OUR PRODUCTS?  
WE WILL GLADLY  
SUBMIT YOU  
AN OFFER.

(if available)

Quantity

WWW.STEMA-MEDIZINTECHNIK.DE

## 7. INDEX

|               |                                                              |               |                                         |               |                                         |
|---------------|--------------------------------------------------------------|---------------|-----------------------------------------|---------------|-----------------------------------------|
| AA-002-001    | <a href="#">8</a>                                            | AC-002-004    | <a href="#">33</a>                      | AC-004-204    | <a href="#">35</a>                      |
| AA-002-002    | <a href="#">8</a>                                            | AC-002-004-R  | <a href="#">33</a>                      | AC-004-204-OT | <a href="#">38</a>                      |
| AA-002-003    | <a href="#">8</a> , <a href="#">59</a>                       | AC-002-005    | <a href="#">33</a>                      | AC-004-204-T  | <a href="#">38</a>                      |
| AA-002-011    | <a href="#">8</a>                                            | AC-002-005-R  | <a href="#">33</a>                      | AC-004-205    | <a href="#">35</a>                      |
| AA-002-012    | <a href="#">8</a> , <a href="#">59</a>                       | AC-002-201    | <a href="#">33</a>                      | AC-004-205-OT | <a href="#">38</a>                      |
| AA-002-013    | <a href="#">8</a> , <a href="#">59</a>                       | AC-002-202    | <a href="#">33</a>                      | AC-004-205-T  | <a href="#">38</a>                      |
| AA-002-013    | <a href="#">77</a> , <a href="#">78</a>                      | AC-002-203    | <a href="#">33</a>                      | AC-005-004    | <a href="#">35</a>                      |
| AA-002-014    | <a href="#">8</a>                                            | AC-002-204    | <a href="#">33</a>                      | AC-005-004-R  | <a href="#">35</a>                      |
| AA-002-023    | <a href="#">59</a>                                           | AC-002-205    | <a href="#">33</a>                      | AC-005-204    | <a href="#">35</a>                      |
| AA-003-001    | <a href="#">17</a>                                           | AC-002-502    | <a href="#">36</a>                      | AC-006-001    | <a href="#">33</a>                      |
| AA-003-002    | <a href="#">17</a> , <a href="#">65</a> , <a href="#">70</a> | AC-002-512    | <a href="#">36</a>                      | AC-006-001-R  | <a href="#">33</a>                      |
| AA-003-004    | <a href="#">65</a>                                           | AC-003-000    | <a href="#">35</a>                      | AC-006-002    | <a href="#">33</a> , <a href="#">79</a> |
| AA-003-011    | <a href="#">17</a>                                           | AC-003-000-R  | <a href="#">35</a>                      | AC-006-002-R  | <a href="#">33</a>                      |
| AA-003-012    | <a href="#">17</a> , <a href="#">65</a> , <a href="#">70</a> | AC-003-001    | <a href="#">35</a>                      | AC-006-003    | <a href="#">33</a>                      |
| AA-003-012    | <a href="#">79</a> , <a href="#">80</a>                      | AC-003-001-R  | <a href="#">35</a>                      | AC-006-003-R  | <a href="#">33</a>                      |
| AA-003-014    | <a href="#">17</a> , <a href="#">65</a>                      | AC-003-002    | <a href="#">35</a> , <a href="#">79</a> | AC-006-004    | <a href="#">33</a> , <a href="#">78</a> |
| AA-004-001    | <a href="#">26</a>                                           | AC-003-002-R  | <a href="#">35</a>                      | AC-006-004-R  | <a href="#">33</a>                      |
| AA-004-002    | <a href="#">26</a>                                           | AC-003-003    | <a href="#">35</a>                      | AC-006-005    | <a href="#">33</a>                      |
| AA-004-003    | <a href="#">26</a>                                           | AC-003-003-R  | <a href="#">35</a>                      | AC-006-005-R  | <a href="#">33</a>                      |
| AC-001-001    | <a href="#">33</a>                                           | AC-003-004    | <a href="#">35</a> , <a href="#">78</a> | AC-006-201    | <a href="#">33</a>                      |
| AC-001-001-R  | <a href="#">33</a>                                           | AC-003-004-R  | <a href="#">35</a>                      | AC-006-201-OT | <a href="#">38</a>                      |
| AC-001-002    | <a href="#">33</a> , <a href="#">79</a>                      | AC-003-005    | <a href="#">35</a>                      | AC-006-201-T  | <a href="#">38</a>                      |
| AC-001-002-R  | <a href="#">33</a>                                           | AC-003-005-R  | <a href="#">35</a>                      | AC-006-202    | <a href="#">33</a>                      |
| AC-001-003    | <a href="#">33</a>                                           | AC-003-200    | <a href="#">35</a>                      | AC-006-202-OT | <a href="#">38</a>                      |
| AC-001-003-R  | <a href="#">33</a>                                           | AC-003-200-OT | <a href="#">38</a>                      | AC-006-202-T  | <a href="#">38</a>                      |
| AC-001-004    | <a href="#">33</a> , <a href="#">78</a>                      | AC-003-200-T  | <a href="#">38</a>                      | AC-006-203    | <a href="#">33</a>                      |
| AC-001-004-R  | <a href="#">33</a>                                           | AC-003-201    | <a href="#">35</a>                      | AC-006-203-OT | <a href="#">38</a>                      |
| AC-001-005    | <a href="#">33</a>                                           | AC-003-201-OT | <a href="#">38</a>                      | AC-006-203-T  | <a href="#">38</a>                      |
| AC-001-005-R  | <a href="#">33</a>                                           | AC-003-201-T  | <a href="#">38</a>                      | AC-006-204    | <a href="#">33</a>                      |
| AC-001-008    | <a href="#">33</a>                                           | AC-003-202    | <a href="#">35</a>                      | AC-006-204-OT | <a href="#">38</a>                      |
| AC-001-008-R  | <a href="#">33</a>                                           | AC-003-202-OT | <a href="#">38</a>                      | AC-006-204-T  | <a href="#">38</a>                      |
| AC-001-201    | <a href="#">33</a>                                           | AC-003-202-T  | <a href="#">38</a>                      | AC-006-205    | <a href="#">33</a>                      |
| AC-001-201-OT | <a href="#">38</a>                                           | AC-003-203    | <a href="#">35</a>                      | AC-006-205-OT | <a href="#">38</a>                      |
| AC-001-201-T  | <a href="#">38</a>                                           | AC-003-204    | <a href="#">35</a>                      | AC-006-205-T  | <a href="#">38</a>                      |
| AC-001-202    | <a href="#">33</a>                                           | AC-003-204-OT | <a href="#">38</a>                      | AC-006-502    | <a href="#">36</a>                      |
| AC-001-202-OT | <a href="#">38</a>                                           | AC-003-204-T  | <a href="#">38</a>                      | AC-006-504    | <a href="#">36</a>                      |
| AC-001-202-T  | <a href="#">38</a>                                           | AC-003-205    | <a href="#">35</a>                      | AC-006-512    | <a href="#">36</a>                      |
| AC-001-203    | <a href="#">33</a>                                           | AC-003-205-OT | <a href="#">38</a>                      | AC-006-514    | <a href="#">36</a>                      |
| AC-001-203-OT | <a href="#">38</a>                                           | AC-003-205-T  | <a href="#">38</a>                      | AC-007-001    | <a href="#">34</a>                      |
| AC-001-203-T  | <a href="#">38</a>                                           | AC-003-302    | <a href="#">35</a>                      | AC-007-001-R  | <a href="#">34</a>                      |
| AC-001-204    | <a href="#">33</a>                                           | AC-003-304    | <a href="#">35</a>                      | AC-007-002    | <a href="#">34</a>                      |
| AC-001-204-OT | <a href="#">38</a>                                           | AC-003-502    | <a href="#">36</a>                      | AC-007-002-R  | <a href="#">34</a>                      |
| AC-001-204-T  | <a href="#">38</a>                                           | AC-003-504    | <a href="#">36</a>                      | AC-007-201    | <a href="#">34</a>                      |
| AC-001-205    | <a href="#">33</a>                                           | AC-003-512    | <a href="#">36</a>                      | AC-007-201-OT | <a href="#">38</a>                      |
| AC-001-206-OT | <a href="#">38</a>                                           | AC-003-514    | <a href="#">36</a>                      | AC-007-201-T  | <a href="#">38</a>                      |
| AC-001-206-T  | <a href="#">38</a>                                           | AC-003-602    | <a href="#">36</a>                      | AC-007-202    | <a href="#">34</a>                      |
| AC-001-208    | <a href="#">33</a>                                           | AC-003-604    | <a href="#">36</a>                      | AC-007-202-OT | <a href="#">38</a>                      |
| AC-001-502    | <a href="#">36</a>                                           | AC-003-612    | <a href="#">36</a>                      | AC-007-202-T  | <a href="#">38</a>                      |
| AC-001-504    | <a href="#">36</a>                                           | AC-003-612-OT | <a href="#">38</a>                      | AC-008-004    | <a href="#">34</a>                      |
| AC-001-512    | <a href="#">36</a>                                           | AC-003-612-T  | <a href="#">38</a>                      | AC-008-004-R  | <a href="#">34</a>                      |
| AC-001-514    | <a href="#">36</a>                                           | AC-003-614    | <a href="#">36</a>                      | AC-008-005    | <a href="#">34</a>                      |
| AC-002-001    | <a href="#">33</a>                                           | AC-003-614-OT | <a href="#">38</a>                      | AC-008-005-R  | <a href="#">34</a>                      |
| AC-002-001-R  | <a href="#">33</a>                                           | AC-003-614-T  | <a href="#">38</a>                      | AC-008-204    | <a href="#">34</a>                      |
| AC-002-002    | <a href="#">33</a>                                           | AC-004-004    | <a href="#">35</a>                      | AC-008-205    | <a href="#">34</a>                      |
| AC-002-002-R  | <a href="#">33</a>                                           | AC-004-004-R  | <a href="#">35</a>                      | AC-009-001    | <a href="#">34</a>                      |
| AC-002-003    | <a href="#">33</a>                                           | AC-004-005    | <a href="#">35</a>                      | AC-009-001-R  | <a href="#">34</a>                      |
| AC-002-003-R  | <a href="#">33</a>                                           | AC-004-005-R  | <a href="#">35</a>                      | AC-009-002    | <a href="#">34</a>                      |

## 7. INDEX

|               |                    |               |                        |            |                        |
|---------------|--------------------|---------------|------------------------|------------|------------------------|
| AC-009-002-R  | <a href="#">34</a> | AC-020-004-R  | <a href="#">34</a>     | AD-001-013 | <a href="#">40</a>     |
| AC-009-004    | <a href="#">34</a> | AC-020-202    | <a href="#">34</a>     | AD-001-014 | <a href="#">40</a>     |
| AC-009-004-R  | <a href="#">34</a> | AC-020-204    | <a href="#">34</a>     | AD-001-015 | <a href="#">40</a>     |
| AC-009-201    | <a href="#">34</a> | AC-101-201-OT | <a href="#">39</a>     | AD-001-016 | <a href="#">40</a>     |
| AC-009-201-OT | <a href="#">38</a> | AC-101-201-T  | <a href="#">39</a>     | AD-001-017 | <a href="#">40</a>     |
| AC-009-201-T  | <a href="#">38</a> | AC-101-202-OT | <a href="#">39</a>     | AD-001-018 | <a href="#">40</a>     |
| AC-009-202    | <a href="#">34</a> | AC-101-202-T  | <a href="#">39</a>     | AD-002-001 | <a href="#">40</a>     |
| AC-009-202-OT | <a href="#">38</a> | AC-101-203-OT | <a href="#">39</a>     | AD-002-002 | <a href="#">40</a>     |
| AC-009-204    | <a href="#">34</a> | AC-101-203-T  | <a href="#">39</a>     | AD-002-003 | <a href="#">40, 79</a> |
| AC-010-004    | <a href="#">33</a> | AC-101-204-OT | <a href="#">39</a>     | AD-002-005 | <a href="#">40, 78</a> |
| AC-010-204    | <a href="#">33</a> | AC-101-204-T  | <a href="#">39</a>     | AD-002-011 | <a href="#">40</a>     |
| AC-010-502    | <a href="#">36</a> | AC-101-206-OT | <a href="#">39</a>     | AD-002-012 | <a href="#">40</a>     |
| AC-010-504    | <a href="#">36</a> | AC-101-206-T  | <a href="#">39</a>     | AD-002-013 | <a href="#">40</a>     |
| AC-010-512    | <a href="#">36</a> | AC-103-200-OT | <a href="#">39</a>     | AD-002-015 | <a href="#">40</a>     |
| AC-010-514    | <a href="#">36</a> | AC-103-200-T  | <a href="#">39</a>     | AD-003-003 | <a href="#">40</a>     |
| AC-011-004    | <a href="#">33</a> | AC-103-201-OT | <a href="#">39</a>     | AD-003-005 | <a href="#">40</a>     |
| AC-011-204    | <a href="#">33</a> | AC-103-201-T  | <a href="#">39</a>     | AD-003-007 | <a href="#">40</a>     |
| AC-012-001    | <a href="#">34</a> | AC-103-202-OT | <a href="#">39</a>     | AD-003-008 | <a href="#">40</a>     |
| AC-012-002    | <a href="#">34</a> | AC-103-202-T  | <a href="#">39</a>     | AD-003-013 | <a href="#">40</a>     |
| AC-012-004    | <a href="#">34</a> | AC-103-204-OT | <a href="#">39</a>     | AD-003-015 | <a href="#">40</a>     |
| AC-014-502    | <a href="#">36</a> | AC-103-204-T  | <a href="#">39</a>     | AD-003-017 | <a href="#">40</a>     |
| AC-014-504    | <a href="#">36</a> | AC-103-205-OT | <a href="#">39</a>     | AD-003-018 | <a href="#">40</a>     |
| AC-014-512    | <a href="#">36</a> | AC-103-205-T  | <a href="#">39</a>     | AD-004-001 | <a href="#">43</a>     |
| AC-014-514    | <a href="#">36</a> | AC-103-612-OT | <a href="#">39</a>     | AD-004-002 | <a href="#">43</a>     |
| AC-015-002    | <a href="#">33</a> | AC-103-612-T  | <a href="#">39</a>     | AD-004-003 | <a href="#">43</a>     |
| AC-015-002-R  | <a href="#">33</a> | AC-103-614-OT | <a href="#">39</a>     | AD-004-004 | <a href="#">43</a>     |
| AC-015-003-R  | <a href="#">33</a> | AC-103-614-T  | <a href="#">39</a>     | AD-010-001 | <a href="#">41</a>     |
| AC-015-004    | <a href="#">33</a> | AC-104-204-OT | <a href="#">39</a>     | AD-010-002 | <a href="#">41</a>     |
| AC-015-004-R  | <a href="#">33</a> | AC-104-204-T  | <a href="#">39</a>     | AD-010-003 | <a href="#">41</a>     |
| AC-015-202    | <a href="#">33</a> | AC-104-205-OT | <a href="#">39</a>     | AD-010-004 | <a href="#">41</a>     |
| AC-015-203    | <a href="#">33</a> | AC-104-205-T  | <a href="#">39</a>     | AD-010-005 | <a href="#">41</a>     |
| AC-015-204    | <a href="#">33</a> | AC-106-201-OT | <a href="#">39</a>     | AD-010-006 | <a href="#">41</a>     |
| AC-016-001    | <a href="#">34</a> | AC-106-201-T  | <a href="#">39</a>     | AD-010-007 | <a href="#">41</a>     |
| AC-016-001-R  | <a href="#">34</a> | AC-106-202-OT | <a href="#">39</a>     | AD-010-008 | <a href="#">41</a>     |
| AC-016-002    | <a href="#">34</a> | AC-106-202-T  | <a href="#">39</a>     | AD-010-011 | <a href="#">75</a>     |
| AC-016-002-R  | <a href="#">34</a> | AC-106-203-OT | <a href="#">39</a>     | AD-010-012 | <a href="#">75</a>     |
| AC-016-201    | <a href="#">34</a> | AC-106-203-T  | <a href="#">39</a>     | AD-010-013 | <a href="#">75</a>     |
| AC-016-202    | <a href="#">34</a> | AC-106-204-OT | <a href="#">39</a>     | AD-010-015 | <a href="#">75</a>     |
| AC-017-001    | <a href="#">33</a> | AC-106-204-T  | <a href="#">39</a>     | AD-010-024 | <a href="#">41</a>     |
| AC-017-001-R  | <a href="#">33</a> | AC-106-205-OT | <a href="#">39</a>     | AD-010-025 | <a href="#">42</a>     |
| AC-017-004    | <a href="#">33</a> | AC-106-205-T  | <a href="#">39</a>     | AD-010-026 | <a href="#">42</a>     |
| AC-017-004-R  | <a href="#">33</a> | AC-107-201-OT | <a href="#">39</a>     | AE-001-030 | <a href="#">63</a>     |
| AC-017-201    | <a href="#">33</a> | AC-107-201-T  | <a href="#">39</a>     | AE-001-031 | <a href="#">63</a>     |
| AC-017-204    | <a href="#">33</a> | AC-107-202-OT | <a href="#">39</a>     | AE-001-032 | <a href="#">63</a>     |
| AC-018-001-R  | <a href="#">35</a> | AC-107-202-T  | <a href="#">39</a>     | AE-001-033 | <a href="#">63</a>     |
| AC-018-002-R  | <a href="#">35</a> | AC-109-201-OT | <a href="#">39</a>     | AE-001-036 | <a href="#">63</a>     |
| AC-018-201    | <a href="#">35</a> | AC-109-201-T  | <a href="#">39</a>     | AE-001-037 | <a href="#">63</a>     |
| AC-018-202    | <a href="#">35</a> | AD-001-001    | <a href="#">40</a>     | AE-001-038 | <a href="#">63</a>     |
| AC-019-001    | <a href="#">34</a> | AD-001-002    | <a href="#">40</a>     | AE-001-039 | <a href="#">64</a>     |
| AC-019-001-R  | <a href="#">34</a> | AD-001-003    | <a href="#">40, 79</a> | AE-001-040 | <a href="#">63</a>     |
| AC-019-002    | <a href="#">34</a> | AD-001-004    | <a href="#">40</a>     | AE-001-041 | <a href="#">63</a>     |
| AC-019-002-R  | <a href="#">34</a> | AD-001-005    | <a href="#">40, 78</a> | AE-001-042 | <a href="#">63</a>     |
| AC-019-201    | <a href="#">34</a> | AD-001-006    | <a href="#">40</a>     | AE-001-043 | <a href="#">63</a>     |
| AC-019-202    | <a href="#">34</a> | AD-001-007    | <a href="#">40</a>     | AE-001-044 | <a href="#">64</a>     |
| AC-020-002    | <a href="#">34</a> | AD-001-008    | <a href="#">40</a>     | AE-001-045 | <a href="#">64</a>     |
| AC-020-002-R  | <a href="#">34</a> | AD-001-011    | <a href="#">40</a>     | AE-001-046 | <a href="#">64</a>     |
| AC-020-004    | <a href="#">34</a> | AD-001-012    | <a href="#">40</a>     | AE-001-047 | <a href="#">64</a>     |

## 7. INDEX

|            |                        |            |                    |            |                    |
|------------|------------------------|------------|--------------------|------------|--------------------|
| AE-001-048 | <a href="#">63</a>     | AJ-003-009 | <a href="#">50</a> | AJ-003-074 | <a href="#">50</a> |
| AE-002-031 | <a href="#">64</a>     | AJ-003-010 | <a href="#">50</a> | AJ-003-078 | <a href="#">54</a> |
| AE-003-031 | <a href="#">69</a>     | AJ-003-011 | <a href="#">50</a> | AJ-003-079 | <a href="#">54</a> |
| AE-003-032 | <a href="#">69</a>     | AJ-003-012 | <a href="#">50</a> | AJ-003-081 | <a href="#">52</a> |
| AE-003-033 | <a href="#">69</a>     | AJ-003-013 | <a href="#">50</a> | AJ-003-082 | <a href="#">52</a> |
| AE-003-034 | <a href="#">69</a>     | AJ-003-014 | <a href="#">50</a> | AJ-003-084 | <a href="#">52</a> |
| AE-003-036 | <a href="#">69</a>     | AJ-003-015 | <a href="#">50</a> | AJ-003-086 | <a href="#">52</a> |
| AE-003-040 | <a href="#">69</a>     | AJ-003-016 | <a href="#">50</a> | AJ-003-088 | <a href="#">52</a> |
| AE-003-041 | <a href="#">74, 80</a> | AJ-003-021 | <a href="#">50</a> | AJ-003-090 | <a href="#">52</a> |
| AE-003-043 | <a href="#">74</a>     | AJ-003-022 | <a href="#">50</a> | AJ-003-092 | <a href="#">52</a> |
| AE-003-044 | <a href="#">74, 80</a> | AJ-003-023 | <a href="#">50</a> | AJ-003-094 | <a href="#">52</a> |
| AE-003-045 | <a href="#">74</a>     | AJ-003-024 | <a href="#">50</a> | AJ-003-096 | <a href="#">52</a> |
| AE-003-046 | <a href="#">74, 80</a> | AJ-003-025 | <a href="#">50</a> | AJ-003-098 | <a href="#">52</a> |
| AG-001-011 | <a href="#">62</a>     | AJ-003-026 | <a href="#">50</a> | AJ-003-100 | <a href="#">52</a> |
| AG-001-012 | <a href="#">62</a>     | AJ-003-027 | <a href="#">50</a> | AJ-003-102 | <a href="#">52</a> |
| AG-002-011 | <a href="#">68</a>     | AJ-003-028 | <a href="#">50</a> | AJ-003-104 | <a href="#">52</a> |
| AG-002-012 | <a href="#">68</a>     | AJ-003-029 | <a href="#">50</a> | AJ-003-106 | <a href="#">52</a> |
| AG-002-013 | <a href="#">71</a>     | AJ-003-030 | <a href="#">50</a> | AJ-003-108 | <a href="#">52</a> |
| AG-002-015 | <a href="#">71</a>     | AJ-003-031 | <a href="#">50</a> | AJ-003-110 | <a href="#">52</a> |
| AG-002-015 | <a href="#">80</a>     | AJ-003-032 | <a href="#">50</a> | AJ-003-112 | <a href="#">52</a> |
| AH-001-011 | <a href="#">61</a>     | AJ-003-033 | <a href="#">50</a> | AJ-003-114 | <a href="#">52</a> |
| AH-001-012 | <a href="#">61</a>     | AJ-003-034 | <a href="#">50</a> | AJ-003-116 | <a href="#">52</a> |
| AH-001-013 | <a href="#">61</a>     | AJ-003-035 | <a href="#">50</a> | AJ-003-118 | <a href="#">52</a> |
| AH-001-015 | <a href="#">61</a>     | AJ-003-036 | <a href="#">50</a> | AJ-003-120 | <a href="#">52</a> |
| AH-001-016 | <a href="#">61</a>     | AJ-003-037 | <a href="#">50</a> | AJ-003-121 | <a href="#">52</a> |
| AH-001-017 | <a href="#">61</a>     | AJ-003-038 | <a href="#">50</a> | AJ-003-122 | <a href="#">52</a> |
| AH-001-018 | <a href="#">61</a>     | AJ-003-039 | <a href="#">50</a> | AJ-003-123 | <a href="#">52</a> |
| AH-001-301 | <a href="#">60</a>     | AJ-003-040 | <a href="#">50</a> | AJ-003-124 | <a href="#">52</a> |
| AH-001-302 | <a href="#">60</a>     | AJ-003-041 | <a href="#">50</a> | AJ-003-125 | <a href="#">52</a> |
| AH-001-303 | <a href="#">60</a>     | AJ-003-042 | <a href="#">50</a> | AJ-003-126 | <a href="#">52</a> |
| AH-001-304 | <a href="#">60</a>     | AJ-003-043 | <a href="#">50</a> | AJ-003-127 | <a href="#">52</a> |
| AH-001-305 | <a href="#">60</a>     | AJ-003-044 | <a href="#">50</a> | AJ-003-128 | <a href="#">52</a> |
| AH-001-306 | <a href="#">60</a>     | AJ-003-045 | <a href="#">50</a> | AJ-003-151 | <a href="#">51</a> |
| AH-001-311 | <a href="#">60</a>     | AJ-003-046 | <a href="#">50</a> | AJ-003-152 | <a href="#">51</a> |
| AH-001-312 | <a href="#">60</a>     | AJ-003-050 | <a href="#">50</a> | AJ-003-153 | <a href="#">51</a> |
| AH-002-011 | <a href="#">67</a>     | AJ-003-051 | <a href="#">50</a> | AJ-003-154 | <a href="#">51</a> |
| AH-002-013 | <a href="#">67</a>     | AJ-003-052 | <a href="#">50</a> | AJ-003-155 | <a href="#">51</a> |
| AH-002-015 | <a href="#">73</a>     | AJ-003-053 | <a href="#">50</a> | AJ-003-156 | <a href="#">51</a> |
| AH-002-015 | <a href="#">80</a>     | AJ-003-054 | <a href="#">50</a> | AJ-003-157 | <a href="#">51</a> |
| AH-002-016 | <a href="#">73</a>     | AJ-003-055 | <a href="#">50</a> | AJ-003-158 | <a href="#">51</a> |
| AH-002-211 | <a href="#">72</a>     | AJ-003-056 | <a href="#">50</a> | AJ-003-170 | <a href="#">51</a> |
| AH-002-211 | <a href="#">80</a>     | AJ-003-057 | <a href="#">50</a> | AJ-003-171 | <a href="#">51</a> |
| AH-002-213 | <a href="#">72</a>     | AJ-003-058 | <a href="#">50</a> | AJ-003-172 | <a href="#">51</a> |
| AH-002-213 | <a href="#">80</a>     | AJ-003-059 | <a href="#">50</a> | AJ-003-173 | <a href="#">51</a> |
| AH-002-216 | <a href="#">72</a>     | AJ-003-060 | <a href="#">50</a> | AJ-003-174 | <a href="#">51</a> |
| AH-002-301 | <a href="#">66</a>     | AJ-003-061 | <a href="#">50</a> | AJ-003-175 | <a href="#">51</a> |
| AH-002-301 | <a href="#">79</a>     | AJ-003-062 | <a href="#">50</a> | AJ-003-176 | <a href="#">51</a> |
| AH-002-302 | <a href="#">66</a>     | AJ-003-063 | <a href="#">50</a> | AJ-003-177 | <a href="#">51</a> |
| AH-003-301 | <a href="#">78</a>     | AJ-003-064 | <a href="#">50</a> | AJ-003-178 | <a href="#">51</a> |
| AJ-003-001 | <a href="#">50</a>     | AJ-003-065 | <a href="#">50</a> | AJ-003-179 | <a href="#">51</a> |
| AJ-003-002 | <a href="#">50</a>     | AJ-003-066 | <a href="#">50</a> | AJ-003-180 | <a href="#">51</a> |
| AJ-003-003 | <a href="#">50</a>     | AJ-003-067 | <a href="#">50</a> | AJ-003-181 | <a href="#">51</a> |
| AJ-003-004 | <a href="#">50</a>     | AJ-003-068 | <a href="#">50</a> | AJ-003-182 | <a href="#">51</a> |
| AJ-003-005 | <a href="#">50</a>     | AJ-003-069 | <a href="#">50</a> | AJ-003-183 | <a href="#">51</a> |
| AJ-003-006 | <a href="#">50</a>     | AJ-003-070 | <a href="#">50</a> | AJ-003-184 | <a href="#">51</a> |
| AJ-003-007 | <a href="#">50</a>     | AJ-003-071 | <a href="#">50</a> | AJ-003-185 | <a href="#">51</a> |
| AJ-003-008 | <a href="#">50</a>     | AJ-003-072 | <a href="#">50</a> | AJ-003-186 | <a href="#">51</a> |

## 7. INDEX

|            |                            |            |                        |              |                        |
|------------|----------------------------|------------|------------------------|--------------|------------------------|
| AJ-003-187 | <a href="#">51</a>         | AM-002-020 | <a href="#">48</a>     | AN-003-002   | <a href="#">9</a>      |
| AJ-003-188 | <a href="#">51</a>         | AM-002-021 | <a href="#">48</a>     | AN-003-003   | <a href="#">9</a>      |
| AJ-003-189 | <a href="#">51</a>         | AM-002-022 | <a href="#">48</a>     | AN-003-101   | <a href="#">10</a>     |
| AJ-003-190 | <a href="#">51</a>         | AM-002-023 | <a href="#">48</a>     | AN-003-102   | <a href="#">10, 77</a> |
| AJ-003-191 | <a href="#">51</a>         | AM-002-024 | <a href="#">48</a>     | AN-003-103   | <a href="#">10, 77</a> |
| AJ-003-192 | <a href="#">51</a>         | AM-002-025 | <a href="#">48</a>     | AN-003-104   | <a href="#">10, 77</a> |
| AJ-003-193 | <a href="#">51</a>         | AM-002-026 | <a href="#">48</a>     | AN-003-201   | <a href="#">11</a>     |
| AJ-003-194 | <a href="#">51</a>         | AM-002-027 | <a href="#">48</a>     | AN-003-202   | <a href="#">11</a>     |
| AJ-003-195 | <a href="#">51</a>         | AM-002-028 | <a href="#">48</a>     | AN-003-301   | <a href="#">12</a>     |
| AJ-003-196 | <a href="#">51</a>         | AM-002-029 | <a href="#">48</a>     | AN-003-302   | <a href="#">12, 78</a> |
| AJ-003-197 | <a href="#">51</a>         | AM-002-030 | <a href="#">48</a>     | AN-003-303   | <a href="#">12</a>     |
| AJ-003-198 | <a href="#">51</a>         | AM-002-031 | <a href="#">48</a>     | AN-004-401   | <a href="#">14</a>     |
| AJ-003-199 | <a href="#">51</a>         | AM-002-032 | <a href="#">48</a>     | AN-004-402   | <a href="#">14</a>     |
| AJ-003-200 | <a href="#">51</a>         | AM-002-040 | <a href="#">54</a>     | AN-004-403   | <a href="#">14</a>     |
| AJ-003-201 | <a href="#">51</a>         | AM-002-041 | <a href="#">54</a>     | AN-005-301   | <a href="#">13</a>     |
| AJ-003-202 | <a href="#">51</a>         | AM-004-001 | <a href="#">45</a>     | AN-005-302   | <a href="#">13</a>     |
| AJ-003-203 | <a href="#">51</a>         | AM-004-002 | <a href="#">45</a>     | AN-005-303   | <a href="#">13</a>     |
| AJ-005-001 | <a href="#">49</a>         | AM-004-004 | <a href="#">45</a>     | AN-005-304   | <a href="#">13</a>     |
| AJ-005-002 | <a href="#">49</a>         | AM-004-006 | <a href="#">45</a>     |              |                        |
| AJ-006-001 | <a href="#">53</a>         | AN-001-001 | <a href="#">27</a>     |              |                        |
| AJ-006-002 | <a href="#">53</a>         | AN-001-002 | <a href="#">27</a>     | CI-001-030   | <a href="#">42</a>     |
| AJ-006-003 | <a href="#">53</a>         | AN-001-003 | <a href="#">27</a>     | CI-001-040   | <a href="#">42</a>     |
| AJ-006-004 | <a href="#">53</a>         | AN-001-101 | <a href="#">28</a>     |              |                        |
| AJ-006-005 | <a href="#">53</a>         | AN-001-102 | <a href="#">28</a>     |              |                        |
| AJ-006-006 | <a href="#">53</a>         | AN-001-103 | <a href="#">28</a>     |              |                        |
| AJ-006-007 | <a href="#">53</a>         | AN-001-104 | <a href="#">28</a>     | EK-010-010   | <a href="#">44</a>     |
| AJ-006-008 | <a href="#">53</a>         | AN-001-301 | <a href="#">29</a>     | EK-010-011   | <a href="#">44</a>     |
| AJ-006-012 | <a href="#">54</a>         | AN-001-302 | <a href="#">29</a>     | EK-010-012   | <a href="#">44</a>     |
| AJ-010-006 | <a href="#">64, 69, 74</a> | AN-001-310 | <a href="#">30</a>     | EK-010-013   | <a href="#">44</a>     |
| AM-001-001 | <a href="#">45</a>         | AN-001-311 | <a href="#">30</a>     | EK-010-014   | <a href="#">44</a>     |
| AM-001-002 | <a href="#">45</a>         | AN-002-001 | <a href="#">18</a>     | EK-010-015   | <a href="#">44</a>     |
| AM-001-003 | <a href="#">45</a>         | AN-002-002 | <a href="#">18</a>     | EK-010-016   | <a href="#">44</a>     |
| AM-001-004 | <a href="#">45</a>         | AN-002-003 | <a href="#">18</a>     | EK-010-017   | <a href="#">44</a>     |
| AM-002-001 | <a href="#">46</a>         | AN-002-020 | <a href="#">48</a>     |              |                        |
| AM-002-002 | <a href="#">46</a>         | AN-002-101 | <a href="#">19</a>     |              |                        |
| AM-002-003 | <a href="#">46</a>         | AN-002-102 | <a href="#">19</a>     | FA-001-000   | <a href="#">56</a>     |
| AM-002-004 | <a href="#">46</a>         | AN-002-103 | <a href="#">19</a>     | FA-001-001   | <a href="#">56</a>     |
| AM-002-005 | <a href="#">46</a>         | AN-002-104 | <a href="#">19</a>     | FA-001-001-G | <a href="#">56</a>     |
| AM-002-006 | <a href="#">46</a>         | AN-002-201 | <a href="#">20</a>     | FA-001-001-R | <a href="#">56</a>     |
| AM-002-007 | <a href="#">46</a>         | AN-002-202 | <a href="#">20</a>     | FA-002-001   | <a href="#">56</a>     |
| AM-002-008 | <a href="#">46</a>         | AN-002-301 | <a href="#">21</a>     | FA-002-001-G | <a href="#">56</a>     |
| AM-002-010 | <a href="#">47</a>         | AN-002-302 | <a href="#">21, 79</a> | FA-002-001-R | <a href="#">56</a>     |
| AM-002-011 | <a href="#">47</a>         | AN-002-307 | <a href="#">21</a>     | FA-010-001   | <a href="#">55</a>     |
| AM-002-012 | <a href="#">47</a>         | AN-002-308 | <a href="#">21</a>     | FA-010-002   | <a href="#">55</a>     |
| AM-002-013 | <a href="#">47</a>         | AN-002-310 | <a href="#">23</a>     | FA-010-004   | <a href="#">55</a>     |
| AM-002-014 | <a href="#">47</a>         | AN-002-311 | <a href="#">23</a>     | FA-010-005   | <a href="#">55</a>     |
| AM-002-015 | <a href="#">47</a>         | AN-002-312 | <a href="#">22</a>     | FA-010-006   | <a href="#">55</a>     |
| AM-002-016 | <a href="#">47</a>         | AN-002-313 | <a href="#">22</a>     |              |                        |
| AM-002-018 | <a href="#">47</a>         | AN-003-001 | <a href="#">9</a>      |              |                        |

WE COVER THE FOLLOWING AREAS IN SURGERY WITH HIGH QUALITY PRODUCTS

ARTHROSCOPY

BRONCHOSCOPY

GASTROENTEROLOGY

GYNAECOLOGY

ENDOSCOPY

ENT

LAMINECTOMY

LAPAROSCOPY

MICROMOTOR SYSTEM

OMS

PROCTOSCOPY

THORACOSCOPY

UROLOGY

SPINE

CONTAINER SYSTEMS

STERILIZATION



STEMA MEDIZINTECHNIK GMBH | TAKE-OFF GEWERBEPARK 74 | 78579 NEUHAUSEN OB ECK  
INFO@STEMA-MEDIZINTECHNIK.DE | WWW.STEMA-MEDIZINTECHNIK.DE | +49 7461 910 55 09

10/24 | PRINTED IN GERMANY | COPYRIGHT BY STEMA MEDIZINTECHNIK GMBH |  
WE RESERVE THE RIGHT TO MAKE ALTERATIONS